

Product datasheet for **RR215546L3V**

Ppef1 (NM_001034935) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ppef1 (NM_001034935) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ppef1
Synonyms:	MGC125115
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001034935
ORF Size:	1920 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR215546).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001034935.1 , NP_001030107.1
RefSeq Size:	2497 bp
RefSeq ORF:	1923 bp
Locus ID:	317498
UniProt ID:	Q3SWT6
Cytogenetics:	Xq14
Gene Summary:	May have a role in the recovery or adaptation response of photoreceptors. May have a role in development (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)