

Product datasheet for **RR215564L3V**

Atp5g1 (Atp5mc1) (NM_017311) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Atp5g1 (Atp5mc1) (NM_017311) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Atp5mc1
Synonyms:	Atp5g1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_017311
ORF Size:	408 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR215564).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017311.1 , NP_059007.1
RefSeq Size:	561 bp
RefSeq ORF:	411 bp
Locus ID:	29754
UniProt ID:	Q06645
Cytogenetics:	10q26
Gene Summary:	isoform of mitochondril ATP synthase subunit 9; increased expression occurs in response to ethanol consumption [RGD, Feb 2006]



[View online »](#)