

## Product datasheet for **RR215586L3V**

### Ugdh (NM\_031325) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ugdh (NM_031325) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ugdh
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_031325
ORF Size:	1479 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR215586).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_031325.1</a> , <a href="#">NP_112615.1</a>
RefSeq Size:	2318 bp
RefSeq ORF:	1482 bp
Locus ID:	83472
UniProt ID:	<a href="#">O70199</a>
Cytogenetics:	14p11
Gene Summary:	Catalyzes the formation of UDP-alpha-D-glucuronate, a constituent of complex glycosaminoglycans (By similarity). Required for the biosynthesis of chondroitin sulfate and heparan sulfate. Required for embryonic development via its role in the biosynthesis of glycosaminoglycans (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)