

## Product datasheet for **RR215617L4V**

### Krt72 (NM\_001008809) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Krt72 (NM_001008809) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Krt72
Synonyms:	Kb35
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001008809
ORF Size:	1659 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR215617).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001008809.1</a> , <a href="#">NP_001008809.1</a>
RefSeq Size:	1662 bp
RefSeq ORF:	1662 bp
Locus ID:	406227
UniProt ID:	<a href="#">Q6IG04</a>
Cytogenetics:	7q36
Gene Summary:	Has a role in hair formation. Specific component of keratin intermediate filaments in the inner root sheath (IRS) of the hair follicle (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)