

Product datasheet for **RR215623**

Nmnat3 (NM_001013224) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Nmnat3 (NM_001013224) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Nmnat3
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR215623 representing NM_001013224
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAGAACCGAATCCCCGTGGTCTCCTGGCCTGTGGTTCCTCAACCCATCACGAATATGCACCTGC
GTTTGTGGAGGTGGCCAGAGACCCTGCACCAGACAGGAAAGTACCAGGTGATCGAGGGCATCATCTC
GCCTGTCAATGACAGTTATGGGAAGAAAGACCTGGTGGCTTCCACCACCGCTGGCCATGGCCCGCTG
GCCCTGCAGACATCTGACTGGATTCGGGTAGACCCTGGGAGAGCGAGCAGGCACAGTGGATGAAACGG
CATCATCACGGGAGCTGCTCAGATCCGTGGCCAGATGGATGGCCAGACCCAGCAAATACCATCAG
CCTCTGCGGCACTGCC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR215623 representing NM_001013224
Red=Cloning site Green=Tags(s)

MKNRIPVLLACGSFNPITNMHLRLEFVARDHLHQTKYQVIEGII SPVND SYGKKDLVASHHRVAMARL
ALQTSDWIRVDPWESEQAQMETSRRGAAQIRGPDGWPRPQQNTISLCGTA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001013224.1</u> , <u>NP_001013242.1</u>
RefSeq Size:	1854 bp
RefSeq ORF:	369 bp
Locus ID:	363118
Cytogenetics:	8q31
MW:	13.7 kDa