

Product datasheet for **RR215629L4V**

Nampt (NM_177928) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nampt (NM_177928) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Nampt
Synonyms:	Pbef; Pbef1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_177928
ORF Size:	1473 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR215629).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_177928.3 , NP_808789.1
RefSeq Size:	2590 bp
RefSeq ORF:	1476 bp
Locus ID:	297508
UniProt ID:	Q80Z29
Cytogenetics:	6q16
Gene Summary:	has similarity to cytokines but is not secreted; localized to the nucleus and cytoplasm, changes in subcellular localization may be associated with the cell cycle [RGD, Feb 2006]



[View online »](#)