

Product datasheet for RR215658

P2ry12 (NM_022800) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	P2ry12 (NM_022800) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	P2ry12
Synonyms:	P2y12
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RR215658 representing NM_022800 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGGTGCCTGGTGCCAACGCCACCTCAGCCAACACCACCTCCATTCTGGGACCAGCACCTGTGCA
GCAGAGACTACAAGATCACCCAGGTTCTTCCCATTGCTCTACACTGTCCTGTTTTTGGTGGGCTCAT
CACGAACAGCTTGGCGATGAGGATTTTCTCCAGATCCGCAGTAAATCGAACTTCATATTTTTCTTAAG
AACACGGTCATCTCTGATCTTCTTATGATCCTAACTTTCTTTCAAATTCAGTGATGCCAACTGG
GAGCTGGGCACCTGAGAACCCTGGTGTGCCAAGTCACTTCAGTCACGTTTTACTTCACAATGTACATCAG
TATCTCGTTCCTCGGATTGATAACCATTGACCGATACCTGAAGACCACCAGACCTTTAAAACCTCCAGC
CCCAGCAATCTTTGGTGCGAAGATTCTTCTGTTGCCATCTGGGCCTTCATGTTCTGCTGTCCTGCTGC
CTAACATGATTCTCACCACAGGAGGCCAAAAGATAAGGACATAACGAAATGTTCTTTTGAAGTCGGA
GTTTGGTCTGGTCTGGCAGGATAGTCAATTACATCTGCCAAGTCATTTCTGGATTAATTTTTAATT
GTTATTGTTTGTACAGCCTCATTACAAAAGAGCTCTATAGGTCCTACGTGAGACGAGGGGTTACGCCA
AAGCTCCCAAGAAAAGGGTGAACATCAAGGTTTTCATCATCATTGCTGTGTTCTTCATTTGCTTCGTTCC
CTTCCACTTTGCACGGATTCCCTACACCCTGAGCCAGACACGGGCCGTCTTTGACTGCAATGCCGAGAAC
ACTCTGTTCTACGTGAAGGAGAGCACCCCTGTGGCTGACGTCCTTGAACGCCTGCTGATCCATTCATCT
ATTTCTTTCTTTGCAAGTCTTTCAGAAATTCCTTGATGAGCATGCTGAGGTGCTCAACATCCGGGGCAA
TAAGAAGAAAGGACAGGAAGGTGGGGACCCAAGCGAGGAGACCCCAATG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR215658 representing NM_022800
Red=Cloning site Green=Tags(s)

MEVPGANATSANTTSIPGTSTLCSRDKITQVLFPLLYTVLFFAGLITNSLAMRIFFQIRSKSNFIIFLK
 NTVISDLLMILTFPFKILSDAKLGAGHLRTLVCQVTSVTFYFTMYISISFLGLITIDRYLKTRPFKTSS
 PSNLLGAKILSVAIWAFMFLLSLPNMILTNRRPKDKDITKCSFLKSEFGLVWHEIVNYICQVIFWINFLI
 VIVCYSLITKELYRSYVRTRGSAKAPKKRVNIKVFIIIIVFFICFVPPHFARIPYTLSQTRAVFDCNAEN
 TLFYVKESTLWLTSLNACLDPFIYFFLCKSFRNSLMSMLRCSTSGANKKKKGQEGGDPSEETPM

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/ja3047_c01.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_022800

ORF Size: 1029 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_022800.1](#), [NP_073637.1](#)

RefSeq Size: 1254 bp

RefSeq ORF: 1032 bp

Locus ID: 64803


UniProt ID: [Q9EPX4](#)

Cytogenetics: 2q26

MW: 39 kDa

Gene Summary: receptor activated by platelet ADP release; mediates adenylate cyclase inhibition and is a target of antithrombotic drugs [RGD, Feb 2006]

Product images:



Circular map for RR215658