

Product datasheet for **RR215658L4V**

P2ry12 (NM_022800) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	P2ry12 (NM_022800) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	P2ry12
Synonyms:	P2y12
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_022800
ORF Size:	1029 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR215658).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_022800.1 , NP_073637.1
RefSeq Size:	1254 bp
RefSeq ORF:	1032 bp
Locus ID:	64803
UniProt ID:	Q9EPX4
Cytogenetics:	2q26
Gene Summary:	receptor activated by platelet ADP release; mediates adenylate cyclase inhibition and is a target of antithrombotic drugs [RGD, Feb 2006]



[View online »](#)