

Product datasheet for **RR215660L4V**

Tmem98 (NM_001007672) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tmem98 (NM_001007672) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Tmem98
Synonyms:	MGC93733
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001007672
ORF Size:	678 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR215660).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001007672.1 , NP_001007673.1
RefSeq Size:	1438 bp
RefSeq ORF:	681 bp
Locus ID:	303356
UniProt ID:	Q6AYS5
Cytogenetics:	10q26



[View online »](#)

Gene Summary:

Functions as a negative regulator of MYRF in oligodendrocyte differentiation and myelination. Interacts with the C-terminal of MYRF inhibiting MYRF self-cleavage and N-fragment nuclear translocation. The secreted form promotes differentiation of T helper 1 cells (Th1).
[UniProtKB/Swiss-Prot Function]