

Product datasheet for **RR215679**

Dlx1 (NM_001100531) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Dlx1 (NM_001100531) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Dlx1
Synonyms: MGC189409
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR215679 representing NM_001100531
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGC**C

ATGACCATGACCACCATGCCAGAAAGTCTCAACAGCCCGGTGTCGGGCAAGGCTGTGTTTATGGAGTTTG
GGCCGCCCAACCAGCAAATGTCTCCTTCTCCCATGTCCCACGGGCACTACTCCATGCACTGTTTACTC
GGCGGGCCATTTCGAGCCCGACGGTGCCTACAGTTCGGCCTCATCCTTCTCCGACCGCTGGGCTACCC
TAGGTCAACTCGGTGAGCAGCCAGCGTCCAGCCCTACATCAGTTCCGTGCAGTCTACCCCGGAGCG
CCAGCCTCGCCAGAGCCGCTGGAGGACCCAGGGCGGACTCCGAGAAGAGCACCGTGGTGAAGGCGG
AGAAGTGCCTTTAACGGCAAAGGGAAAAAGATCCGTAACCCAGGACAATTTATTCCAGTTTGCAGTTG
CAGGCTTTGAACCGGAGGTTCCAACAACTCAGTACCTAGCTCTGCCTGAGAGGGCGGAGCTGGCGGCT
CCTTGGGACTCACACAGACTCAGGTCAAGATATGGTCCAGAACAACGCTCCAAGTTCAAGAAGCTGAT
GAAGCAAGGCGGGCAGCTCTGGAGGGCAGCGCTGGCGAATGGCAGGGCTTGTCTGCCGGCTCCCA
CCGGTACCACCGGCTGGAATCCGAATTCCTCCTCTGGGAAGGGCTCCGGCAGCAGCGCTGGCTCCTACG
TCCCAGCTACACGTCGTGGTATCCCTCAGCGCACCAAGAAGCTATGCAGCAGCCTCAACTGATG

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >RR215679 representing NM_001100531
 Red=Cloning site Green=Tags(s)

MTMTTTPESLNSPVSGKAVFMEFGPPNQMQSPSPMSHGHYSMHCLHSAGHSQPDGAYSSASSFSRPLGYP
 YVNSVSSHASSPYISSVQSYPGSASLAQSRLDPGADSEKSTVVEGGEVRFNGKGGKIRKPRTIYSSLQL
 QALNRRFQQTQYLALPERAELAASLGLTQTQVKIWFQNKRSKFKKLMKQGGAALGALSALANRALSAGSP
 PVPPGWNPNSSSGKSGSSAGSYVPSYTSWYPSAHQEAMQQPQLM


SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-RsrII

Cloning Scheme:



Plasmid Map:



ACCN: NM_001100531
ORF Size: 765 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001100531.1 , NP_001094001.1
RefSeq Size:	2757 bp
RefSeq ORF:	768 bp
Locus ID:	296500
Cytogenetics:	3q23
MW:	27.3 kDa