

Product datasheet for **RR215680L3V**

Yipf4 (NM_001009712) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Yipf4 (NM_001009712) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Yipf4
Synonyms:	RGD1310157
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001009712
ORF Size:	738 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR215680).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001009712.1 , NP_001009712.1
RefSeq Size:	1847 bp
RefSeq ORF:	741 bp
Locus ID:	362699
UniProt ID:	Q5M7T4
Cytogenetics:	6q13
Gene Summary:	Involved in the maintenance of the Golgi structure.[UniProtKB/Swiss-Prot Function]



[View online »](#)