

Product datasheet for **RR215682L4V**

Eya2 (NM_130427) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Eya2 (NM_130427) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Eya2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_130427
ORF Size:	1593 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR215682).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_130427.1 , NP_569111.1
RefSeq Size:	1596 bp
RefSeq ORF:	1596 bp
Locus ID:	156826
Cytogenetics:	3q42
Gene Summary:	member of the eyes absent (EYA) family of proteins; may play a role in eye development; mouse homolog may act as a transcriptional activator [RGD, Feb 2006]



[View online »](#)