

Product datasheet for **RR215693**

Brd9 (NM_001107453) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Brd9 (NM_001107453) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Brd9
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RR215693 representing NM_001107453
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGGCAAGAAGCACAAGAAGCACAAGGCAGAATGGCGCTCGTCGTACGAAGATTATACAGACACACCAC
 TGGAGAAGCCTCTGAAGCTGGTGTCAAGGTGGGAGGAAGTGAAGTGACAGAGCTCTCAGGATCTGGTCA
 CGACTCCAGCTACTATGATGACCGGTGACACCATGAGCGGAAAGACATAGGGAAAAAGAAAGAAAAAG
 AAGAAGAAGTCAGAGAAGGAAAGCACCTTGATGAGGAGGAGAGGAGGAAGGAAGGAAGAAAGAAAC
 GGAAGCGGGAGAAGAACATTGTGACTCAGAGGGGGAGGCTGATGCCCTTGACCTGGGAAGAAGGTGGA
 GGTGGAACACCCCGAGATCGACCAGCTGAGAATGAGAGCACACCCATCCAGAGGCTTCTGGAGCACTTC
 CTCGCCAGCTACAGAGAAAAGATCCTCATGGATTTTTTGCTTTCTGTTACGGATGCAATTGCTCCTG
 GGTATTCAATGATAATAAAACATCCCATGGACTTTGGCACGATGAAAGACAAGATTGTAGCTAATGAATA
 TAAATCAGTCACAGAATTAAGGCAGATTTCAAATTAATGTGTGATAATGCGATGACGTACAATAGACCA
 GACACCGTGACTACAAATTAGCCAAGAAGATCCTGCACGCGGGCTTTAAGATGATGAGCAAAACAGGCAG
 CTCTCTTGGCAGTGAAGACCCCGCAGTTGAGGAACCTGTTCTGAGGTTGTCCAGTGCAAGTAGAAAAC
 TACCAAGAAATCCAAAAGCCGAGTAGAGAAGTCATCAGCTGCATGTTTGAGCCTGAAGGGAATGCCTGC
 AGCCTGACAGACAGCACAGCAGAGGAGCATGTGCTAGCCCTGGTAGAGCACGCAGCTGATGAGGCTCGGG
 ACAGGATTAACCGGTTTCTCCAGGTGGCAAGATGGGGTACCTGAAGAAGCTTGAGATGGAAGTCTACT
 CTACAGTGTGGTGAATGCAACTGAGCCTGATGCTGATGAGGAGGAGACACACCCTGTGGACCTGAGTTCA
 CTGTCCAGCAAGTTGCTCCAGGTTTACAACATTGGGTTTCAAAGATGAAAGAAGAAATAAAGTCACAT
 TCCTCTCCAGTGCCAGCACTGCATTTCAATGCAGAACAACTCTGTGTTTGGGGACTTGAAGTCGGATGA
 GATGGAGCTTCTGTATTCTGCCTATGGAGATGAGACTGGTGTGCAGTGCGCACTGAGCCTGCAGGAATTC
 GTGAAGGATGCTGGAAGCTACAGCAAGAAAATGGTAGATGACCTCCTGGACCAAAATCACAGGTGGCGATC
 ACTCAAGGATGATCTTCCAGCTGAAGCAGAGGAGGAGCATCCCCATGAGACCTGCAGATGAGATGAAGGT
 TGGGGACCCACTGGGAGATAGTGGTGGCCCTGTTCTGGACTTCATGTCAATGAAACAGTATCCTGATGTC
 TCCCTGGATGTGCCATGCTCAGCTCTCTTGGGAAAGTGAAGAAGGAGCTGGACCATGAGGATAGCCACT
 TGAACCTGGACGAGACAGCCAGGCTCCTGCAGGACTTACATGAAGCACAGCAGAGCGAGGAGGCTCCCG
 GCCATCTCCAACCTTAGCTCTCTGCTCCACTGCCTCTGAGAGGGAGCAGCATCTCCAGGAAGTCTTCT
 CGCCTTAGTGTTGGGGAACAGCCGGATGTACCCATGACCCTTACGAATTCCTTCAATCTCCAGAACCTG
 CAGCTCCTGCCAAGAAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RR215693 representing NM_001107453
 Red=Cloning site Green=Tags(s)


MGKHKHKHKAEWSSYEDYDTPLEKPLKLVKVGSEVTELSGSGHDSSYYDDRSHERERHREKHKHK
 KKKSEKEKHLDEEERRKRKEEKKRKREKEHCDSEGEADAFDPGKKVEVEPPDPRAENESTPIQRLLEHF
 LRQLQRKDPHGFFAFPVTDIAIPGYSMIHKHPMDFGTMKDKIVANEYKSVTEFKADFKLCDNAMTYNRP
 DTVYYKLAKKILHAGFKMMSKQAALLGSEDPAVEEPVPEVVPVQVETTKSKKPSREVISCMFEPEGNAC
 SLTDSTAEHVLALVEHADEARDRINRFLPGGKMGYLLKLDGSLLYSVVNATEPDADEEETHPVDLSS
 LSSKLLPGFTTLGFKDERRNKVTFLLSSASTALSMQNNVFGDLKSDMELLYSAYGDETGVQCALSLQEF
 VKDAGSYSKKMVDLLDQITGGDHSRMIFQLKQRRSIPMRPADEMKGVDPLGDSGGPVLDFFMSMKQYPDV
 SLDVSMSSSLGKVKKELDHEDSHLNLDEARLLQDLHEAQAERGGSRPSSNLSSLSTASEREQHPGSPS
 RLSVGEQPDVTHDPYEFLQSPEPAAPAKN

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



ACCN: NM_001107453

ORF Size: 1767 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001107453.1](#), [NP_001100923.1](#)

RefSeq Size: 2363 bp


RefSeq ORF: 1770 bp

Locus ID: 308067

Cytogenetics: 1p11

MW: 66.1 kDa

Product images:



Circular map for RR215693