

Product datasheet for **RR215697L4V**

Pkm (NM_053297) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pkm (NM_053297) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Pkm
Synonyms:	PK; Pk3; Pkm2; PKM12
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_053297
ORF Size:	1593 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR215697).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_053297.2 , NP_445749.1
RefSeq Size:	1971 bp
RefSeq ORF:	1596 bp
Locus ID:	25630
UniProt ID:	P11980
Cytogenetics:	8q24
Gene Summary:	M1 and M2 isozymes of pyruvate kinase; catalyzes the conversion of ATP and pyruvate to ADP and phosphoenolpyruvate in glycolysis [RGD, Feb 2006]



[View online »](#)