

Product datasheet for **RR215805**

Ncmap (NM_001253918) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ncmap (NM_001253918) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Ncmap
Synonyms: Mp11
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RR215805 representing NM_001253918
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**C

ATGACCACAGCCACCACTCTGGGGGACGCCGTCTTCTCGCTGAACATGACCAGGGGAGAAGATATCCTGT
 ATAAAAGTTCCGAGCCATCGTGGCTGCCATCGTGGTGGTTGTCGTCATCATTGTACCTTGGTGTGAT
 CTTGCTGAAGATGTACAACAGGAGAATGAGGACCAGGCGGGAGCTGGAGCCCAAGAGCCCAAGCCACCT
 GTCCCTCCTGCCCTGGACCCAAACAGCAATGGCAGCCAGCAGCCTGCCGCTGTGACCTCCGACCCTGCCG
 ATGTCCCTGTGGAGACTCGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR215805 representing NM_001253918
 Red=Cloning site Green=Tags(s)

MTTATTLGDAVFSNMTRGEDILYKSSGAIVAAIVVVVVIIVTLVLILLKMYNRRMRTRRELEPKSPKPP
 VPPALDPNSNGSQPAAVTSDPADVPVETR

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI


[View online »](#)

Cloning Scheme:



ACCN: NM_001253918

ORF Size: 300 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


RefSeq: [NM_001253918.1](#), [NP_001240847.1](#)

RefSeq Size: 510 bp

RefSeq ORF: 303 bp

Locus ID: 689826
 UniProt ID: [F1M2Z5](#)
 Cytogenetics: 5q36
 MW: 10.8 kDa
 Gene Summary: Plays a role in myelin formation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RR215805