

## Product datasheet for **RR215905**

### Magohb (NM\_001191105) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Magohb (NM\_001191105) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Magohb  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR215905 representing NM\_001191105  
**Red**=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCTTTGGCTAGTGATTTCTACCTCGGGTACTACGTCGGGCACAAGGGCAAGTTTGGGCACGAGTTTC  
TGGAGTTCGAGTTTCGGCCGACGGGAAGCTTAGATATGCCAACACAGCAATTATAAAAATGATGTGAT  
GATCAGAAAAGAGGCTTATGTCCACAAGAGTGAATGGAAGAAGGATCATTGACGACAGTGAA  
ATTACAAAAGGAGGACGACGCTCTGTGGCCGCTCCTGACAGGGTTGGCCGACGAACTTAAATCGTAA  
TCGGGGATGAACACATTTCTTTACCACATCAAAAATAGGTTCTCTCATTGATGTAATCAGTCAAAGA  
TCCCGAAGGCCTCCGTGTATTTACTATTTGGTGAAGACCTGAAATGTTTAGTGTTCAGCCTCATTGGT  
CTACACTCAAGATTAAGCCGATT

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RR215905 representing NM\_001191105  
**Red**=Cloning site **Green**=Tags(s)

MALASDFYLRYVGHKGFHEFLEFEFRPDGKLRVANNVYKNDVMIRKEAYVHKSVMEEELKRIIDSE  
ITKEDDALWPPDRVGRQLEIIVIGDEHISFTTSKIGSLIDVNQSKDPEGLRVFYLVQDLKCLVSLIG  
LHFKIKPI

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001191105

ORF Size: 444 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001191105.1</a></u> , <u><a href="#">NP_001178034.1</a></u>
<b>RefSeq Size:</b>	630 bp
<b>RefSeq ORF:</b>	447 bp
<b>Locus ID:</b>	690303
<b>Cytogenetics:</b>	4q42
<b>MW:</b>	17.3 kDa