

Product datasheet for **RR215915**

Eif5a2 (NM_001100697) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Eif5a2 (NM_001100697) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Eif5a2
Synonyms: eIF5AII
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR215915 representing NM_001100697
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCAGACGAAATTGATTTCACTACTGGAGATGCAGGGCTTCCAGCACTTATCCTATGCAGTGCTCAG
CCCTGCGCAAAAACGGCTTCGTGGTCTCAAAGGGCGTCCATGCAAAATCGTGGAGATGTCAACTCCAA
GACTGGGAAGCATGGTCATGCCAAGGTACACCTTGTGGAAATTGACATTTCACTGGCAAAAATATGAA
GATATTTGTCCTTCCACTCACAACATGGATGTTCCAAACATTAAGAGAAATGATTATCAACTCATATGCA
TTCAGGATGGTTACCTTTCCCTGCTGACAGAACTGGCGAAGTCCGTGAGGATCTCAAGCTGCCAGAAGG
TGAAGTGGCAAAAGAAATCGAAGGGAAGTACAATGCAGGTGAAGATGTGCAGGTGTCCGTATGTGCCCC
ATGAGTGAAGAATACGCTGTAGCCATAAAGCCCTGCAAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR215915 representing NM_001100697
Red=Cloning site Green=Tags(s)

MADEIDFTTGDAGASSTYPMQCSALRKNFVVLKGRPCKIVEMSTSKTGKHGAKVHLVGIDIFTGKKYE
DICPSTHNMDVPNIKRNQYQLICIQDGYLSLLTETGEVREDLKLPEGELGKEIEGKYNAGEDVQVSVMCA
MSEYVAIAIKPK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001100697

ORF Size: 459 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001100697.2, NP_001094167.1</u>
RefSeq Size:	4871 bp
RefSeq ORF:	462 bp
Locus ID:	310261
Cytogenetics:	2q24
MW:	16.8 kDa