

Product datasheet for **RR215956**

Ascl3 (NM_001271234) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ascl3 (NM_001271234) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Ascl3
Synonyms: Sgn1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR215956 representing NM_001271234
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGACACCAGAAGTTATCCCAGCCCTCCAGACAGGCTCTGGTCTTCGACAGTCTGCCCACTTGCCAC
TGTCAGGCCCTTCTACCTGGACCCCATGGTACTGTCCACCTATGCCCGACGCCCGGTGCCGCCCC
TTACACGGACGAGCTGCCTCTGCTGCCCTTTCCAGTGACACCCTGATCATGGACAATTACGGGGACCCA
TACCCCTTCCCTTCCCTATGCCTTACACCAACTACAGGGCGTGTGACTATACCTATGGCCGGCCTTCA
TCCGCAAGAGGAACGAGCGAGAGAGGCAGCGAGTCAAGTGTGAACGAAGGCTATGCCCGGCTGCGACG
GCACCTGCCAGAGGACTACCTGGAGAAACGACTCAGCAAAGTGGAGACCCTCAGAGCTGCCATCAAGTAT
ATCAGTACCTGCAGTCTCTTGTACCCGGATGAGACTGAGACCAAGAAAAACCCTCGAACAACCAACT
GTGGCTCCCTGGACCCGGCTCTGAGAGTCATT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR215956 representing NM_001271234
Red=Cloning site Green=Tags(s)

MDTRSYPSPDRLSVFAESAHLPLSRPFYLDPMVTVHLCPDAPVPAPYTDLPLLPFSSDTLIMDNYGDP
YPFPPMPYTYNRRCDYTYGPAFIRKRNERRQRVKCVNEGYARLRRHLPEDYLEKRLSKVETLRAAIKY
ISYLSLLYPDETETKKNPRTTNCGLDPALRVI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites: SgfI-MluI



Cloning Scheme:



Plasmid Map:



ACCN: NM_001271234

ORF Size: 522 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001271234.1</u> , <u>NP_001258163.1</u>
RefSeq Size:	689 bp
RefSeq ORF:	525 bp
Locus ID:	246301
Cytogenetics:	1q33
MW:	20.3 kDa
Gene Summary:	mouse homolog represses transcription and may play a role in salivary gland development [RGD, Feb 2006]