

Product datasheet for **RR215967**

Obp3 (NM_147215) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Obp3 (NM_147215) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Obp3
Synonyms:	MGC108576
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RR215967 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAAGCTGTTGCTGCTGCTGTGTCTGGCCTGACCCTGGTCTGTGGCCATGCAGAAGAAGCTAGTT
 TCGAGAGAGGGAACCTCGATGTGGACAAGCTCAATGGGATTGGTTTCTATTGTCGTGGCCTCTGATAA
 AAGAGAAAAGATAGAAGAGAACGGCAGCATGAGAGTTTTGTGCAGCATATCGATGCTTGGAGAATTCC
 TTAGGCTTCACGTTCCGTATTAAGGAAAATGGAGTGTGCACAGAATTTCTTTGGTTGCCGACAAAACAG
 CAAAGGATGGCGAATATTTTGTGAGTATGACGGAGAGAATACATTTACTATACTTAAGACAGACTATGA
 CAATTATGTCATGTTTCATCTCGTTAATGTCAACAACGGGAAACCTCCAGCTGATGGAGCTCTATGGC
 AGAACAAAGGATCTGAGTTCAGACATCAAGGAAAAGTTTGCAAACTATGTGTGGCACATGGAATCACTA
 GGGACAATATCATTGACCTAACCAAGACTGATCGCTGTCTCCAGGCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:	>RR215967 protein sequence Red=Cloning site Green=Tags(s)
-------------------	--

MKLLLLLLCLGLTLVCGHAEASFERGNLDVDKLNQDWFISIVVASDKREKIEENGSMRVFVQHIDVLENS
 LGFTFRIKENGVCTEFSLVADKTAKDGEYFVEYDGENFTILKTDYDNYVMFHLVNVNNGETFQLMELYG
 RTKDLSSDIKEKFAKLCVAHGITRDNIIDLTKTDRCLQA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:	SgfI-MluI
--------------------	-----------


[View online »](#)

Cloning Scheme:


ACCN: NM_147215

ORF Size: 537 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_147215.2](#), [NP_671748.2](#)


RefSeq Size: 1006 bp

RefSeq ORF: 540 bp

Locus ID: 259247

Cytogenetics: 5q24
MW: 20.3 kDa
Gene Summary: odorant binding protein and member of the alpha(2u)-globulin family [RGD, Feb 2006]

Product images:



Circular map for RR215967