

## Product datasheet for **RR216019**

### Nt5c (NM\_001271189) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Nt5c (NM\_001271189) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Nt5c  
**Synonyms:** Umph2  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR216019 representing NM\_001271189  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCGGTGAAGCGGCCGGTGC CGCTGCTGGTGGACATGGACGGCGTGCTAGCTGACTTCGAGTCCGGCC  
TTCTGCAGGGTTCCGCCGTCGTTCCCGGGGAACCGCATGTGCCGCTGGAACAGCGCCGGGGCTTCT  
CGCCAGCGAGCAGTACGGAGCCCTGCGGCCGACCTGGCGAAAAAGTGGCCAGTGTCTATGAAGTCCG  
GGCTTTTTCTAACTTGAGCCATCCCTGGAGCCTTGGACGCTTTCGAGAGATGAACGCCATGCAAG  
ACACCGATGTCTTCATCTGCACCACCCACTGCGGAAGTACGACCACTGTGTGGCGGAGAAGTACCGTTG  
GGTGGAGCAGAACCTGGGGCCACAGTTTGTGGAGCGGATTATCCTGACTAGGGACAAGACAGTGGTCATG  
GGAGACCTGCTCATTGATGACAAAGACAACATTCAAGGCCTAGAGGAGACCCCAAGCTGGGAGCACATCT  
TGTTCACTGCTGCCACAATCAGCACCTGGCCCTGCCCCCACTAGGAGACGGCTGCTTCTGGAGTGA  
TGACTGGAGGGGCATCATAGAGAGCAAGCGAGCCAGCCTG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RR216019 representing NM\_001271189  
Red=Cloning site Green=Tags(s)

MAVKRPVVRVLDMDGVLADFESGLLQGFRRRFPGEHPVPLEQRRGFLASEQYALRPDLAEKVASVYEAP  
GFFLNLEPIPGALDALREMNAMQDQDVFICTTPLRKYDHCYGEKYRWVEQNLGPQFVERIILTRDKTVVM  
GDLLIDDKDNIQGLEETPSWEHILFTCCHNQHLALPPTRRRLLSWSDDWRGIIESKRASL

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001271189

ORF Size: 600 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001271189.1](#), [NP\\_001258118.1](#)

**RefSeq Size:** 826 bp

**RefSeq ORF:** 603 bp

**Locus ID:** 100360582

**Cytogenetics:** 10q32.1

**MW:** 22.9 kDa