

Product datasheet for RR216024

Tmem40 (NM_001191573) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tmem40 (NM_001191573) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Tmem40
Synonyms:	RGD1306976
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR216024 representing NM_001191573 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

**ATGGAGGCTTCAGGATCTTCTTCTCAGTCTCCAAACAGCAGTGCAGTCCACAGAGAAACAGAAGATCTAA
ACCACCCAGAGACAGAGTCCACAAGCAACATGGGAAGGCCCGGGAACGCTATGAGAGAGACAAGTCTTC
CTCCTCATCTCCTCCTCCTCTTCTCCTCTTCTCCTCTTCTCCTCTTCTCCTCTCCTCCTCCTCCTCAGT
GCTGACCATGGGAACTGGAGGTGCTGAAGGACGAGCTTCAGCTCTATGGAGGCGCTGCTGGGAAATGG
TGCTTGCTGGGGAGTCAGGACTCCGAAGGAGAGGTTCCGGCTCAGCGGAGGGAGAAGTGGAGGCCTCTCA
GTTAAACAGACTGAACATAAAGAAGGATGATGAGTTTTCCATTTTGTCCTTCTCTGCTTCGCCATTGGA
GCCTTGCTGGTGTGCTACCACTACTATTAGACTGGTTTCATGTCTCTGGGGTGGTCTGCTCACCTTCG
CTTCTCTGGAGACCATCGGCATCTACTTCGGGCTGGTGTACCGGATCCACAGCGTACTTCAGGGCTTCAT
CCCCCTCTTCAGAAGTTCAGGCTGCCAGGTTTCAGGAGAACCAAC**

**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA**

Protein Sequence:	>RR216024 representing NM_001191573 Red=Cloning site Green=Tags(s)
-------------------	---

MEASGSSSQSPNSSAVHRETDLNHPETEFHKQHKGKARERYERDKSSSSSSSSSSSSSSSSSSSSSSSSSS
ADHGELEVLKDELQLYGAAGEMVLAGEGLRRRGSSEGEVEASQLNRLNIKKDDEFHFVLLCFAIG
ALLVCYHYYSDFMSLVGLLTFASLETIGIYFGLVYRIHSLVQGFIPLLQKFRLLPGFRRTN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:	SgfI-MluI
--------------------	-----------




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001191573

ORF Size: 606 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001191573.1, NP_001178502.1

RefSeq Size: 609 bp

RefSeq ORF: 609 bp

Locus ID: 680858

Cytogenetics: 4q42

MW: 22.1 kDa