

Product datasheet for **RR216030**

Stmn4 (NM_001270855) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Stmn4 (NM_001270855) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Stmn4
Synonyms: Lagl; Rb3
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR216030 representing NM_001270855
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGACCCTCGCAGCCTATAAGGAGAAGATGAAGGAACTCCCACTAGTGTCTCTGTTCTGCTCCTGTTTTT
TGCTGATCCCCTGAATAAATCATCCTACAATATGAAGGCTGGTGTGGGAGACAGTGTAGGAGGAAAGG
TCAAAGCCAGCGAAAGGCAGTCTGACTGGAGAGAAAGAAGAGAACAGGCAGATACGGTAGACCTGAAC
TGGTGTGTCATCTCTGATATGGAAGTCATCGAGCTGAATAAGTGTACCTCGGGCCAGTCTTTGAAGTCA
TCCTGAAGCCACCTTCTTTGACGGGGTGCCTGAGTTTAAATGCCTCCCTCCCAAGACGTCGAGACCCATC
GCTAGAAGAGATACAGAAGAAGCTAGAAGCAGCAGAGGAGCGAAGGAAGTACCAGGAAGCTGAGCTCCTA
AAACACCTTGAGAGAAACGAGAGCATGAGCGTGAGGTAATCCAGAAAGCTATCGAGGAAAAACAACCT
TCATCAAGATGGCGAAAGAGAAGCTGGCCAGAAGATGGAGTCCAATAAGGAAAACCGGGAGGCCATCT
GGCTGCCATGTTGGAGCGGCTGCAAGAGAAGGAGCCGCCTGCTGCGCGG

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR216030 representing NM_001270855
Red=Cloning site Green=Tags(s)

MTLAAYKEKMKELPLVSLFCSCFLSDPLNKSSYKYEWCGRQCRRKQSQSRKGSADWRERREQADTVDLN
WCVISDMEVIELNKCTSGQSFEVILKPPSFDGVPFNASLPRRRDPSLEEIQKLEAAEERRKYQEAELL
KHLAEKREHEREVIQKAIENNNFIKMAKEKLAQKMESNKENREAHLAAMLERLQEKEPPAAR

TRTRPLEQKLISEEDLAANDILDYKDDDDK**V**

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001270855

ORF Size: 609 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001270855.1](#), [NP_001257784.1](#)

RefSeq Size: 1414 bp

RefSeq ORF: 612 bp

Locus ID: 79423

Cytogenetics: 15p12

MW: 23.7 kDa

Gene Summary: controls cell proliferation and activities for stathmin [RGD, Feb 2006]