

## Product datasheet for **RR216107**

### **Rbm24 (NM\_001191100) Rat Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Rbm24 (NM\_001191100) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Rbm24  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR216107 representing NM\_001191100  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGCACACCACCCAGAAGGACACGACGTACACCAAGATCTTCGTGGGAGGGCTGCCCTACCACACCACCG  
ACGCCAGCTGCGCAAGTACTTCGAGGTCTTCGGAGACATCGAGGAGGCTGTGGTCATCACCGACCGGCA  
GACGGCAAGTCCCGGGCTACGGATTTGTCACCATGGCTGACCGCGCTGCTGCTGAAAGAGCCTGTAAG  
GATCCTAACCCCATCATCGATGGTAGGAAGGCCAATGTGAATCTGGCGTACTTGGGAGCCAAACCAAGAA  
TCATGCAGCCAGGTTTTGCCTTTGGCGTTCAACAACCTCACCCGCCCCCTATCCAGAGACCTTTTGGGAT  
CCCCGCCACTACGTCTACCCGCAGGCTTTTGTGCAGCCTGGAGTGGTCAATCCCCACGTCCAGCCACG  
GCAGCTGCGGCTTCCACCACACCGTACATAGATTACACTGGAGCCGCTACGCACAGTACTCGGCGGCGG  
CGGCTGCAGCCGCGGGCTGCTGCCTATGACCAGTACCCATATGCAGCATCGCCAGTCCACAGGCTA  
TGTGACCACTGGGGCTACAGCTATGCCGTCCAGCAGCCATCACCGCTGCTGCGCCCGAACGGCTGCT  
GCTGCGGCTGCGGCCGCTGCAGCCGAGCAGCCTTCGGCCAGTATCAGCCTCAACAGCTGCAGACAGACC  
GAATGCAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RR216107 representing NM\_001191100  
**Red=Cloning site Green=Tags(s)**


MHTTQKDTTYTKIFVGGLPYHTTASLRKYFEVFGDIEEAVVITDRQTGKSRGYGFVTMADRAAAERACK  
DPNPIIDGRKANVNLAYLGAKPRIMQPGFAFGVQQLHPALIQRPFIPAHYVYPQAFVQPQGVVIPHVQPT  
AAAASTTPYIDYTGAAAYQYSAAAAAAAAAAAYDQYPYASPAATGYVTTGGYSYAVQPITAAAPGTAA  
AAAAAAAAAAAFGQYQPQLQTDORMQ

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV




Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001191100

ORF Size: 708 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001191100.1, NP_001178029.1</u>
<b>RefSeq Size:</b>	2611 bp
<b>RefSeq ORF:</b>	711 bp
<b>Locus ID:</b>	690139
<b>Cytogenetics:</b>	17p14
<b>MW:</b>	24.9 kDa