

Product datasheet for RR216151

Bdnf (NM_001270637) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Bdnf (NM_001270637) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Bdnf
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR216151 representing NM_001270637
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGATCGCC

ATGACCATCCTTTTCCTTACTATGGTTATTTCACTTTCGGTTGCATGAAGGCTGCGCCCATGAAAGAAG
CAAACGTCACGACGACAAGGCAACTTGGCCTACCCAGCTGTGCGGACCCATGGGACTCTGGAGAGCGTGAA
TGGGCCAGGGCAGGTTTCGAGAGGCTGACGACGACGTCCTGGCTGACACTTTTGAGCAGTGATCGAA
GAGCTGCTGGATGAGGACCAGAAGGTTTCGCCCCAACGAAGAAAACCATAAGGACGCGGACTTGTACTT
CCCGGGTGATGCTCAGCAGTCAAGTGCCTTTGGAGCCTCCTCTGCTCTTTCTGCTGGAGGAATACAAAA
TTACCTGGATGCCGAAACATGTCTATGAGGGTTCGGCGCCACTCCGACCCCGCCCGCGTGGGAGCTG
AGCGTGTGTGACAGTATTAGCGAGTGGGTACAGCGGCAGATAAAAAGACTGCAGTGGACATGTCCGGT
GGACGGTACAGTCTGGAGAAAGTCCCGGTATCAAAAGGCCAACTGAAGCAATATTTCTACGAGACCAA
GTGTAATCCCATGGGTTACACGAAGGAAGGCTGCAGGGGCATAGACAAAAGGCACTGGAATCGCAATGC
CGAACTACCAATCGTATGTTCCGGCCCTTACTATGGATAGCAAAAAGAGAATTGGCTGGCGGTTTATAA
GGATAGACTTCTGTGTATGTACTGACCATTAAGGGGAAGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR216151 representing NM_001270637
Red=Cloning site Green=Tags(s)

MTILFLTMVISYFGCMKAAPMKEANVHGQGNLAYPAVRTHGTLESVNGPRAGSRGLTTTSLADTFEIVIE
ELLDEDQKVRPNEENHKDADLYTSRVMLSSQVPLEPPLFLLEEYKNYLDAAANMSMRVRRHSDPARRGEL
SVCDISSEWVTAADKKTAVDMSGGTVTVLEKVPVSKGQLKQFYETKCNPMGYTKEGCRGIDKRHWNSQC
RTTQSYVRALTMSKRRIGWRFIRIDTSCVCTLTIKGR


TRTRPLEQKLISEEDLAANDILDYKDDDDKV




[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001270637

ORF Size: 747 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001270637.1</u> , <u>NP_001257566.1</u>
RefSeq Size:	3874 bp
RefSeq ORF:	750 bp
Locus ID:	24225
UniProt ID:	<u>P23363</u>
Cytogenetics:	3q34
MW:	28.1 kDa
Gene Summary:	plays a role in the development of hippocampal long term potentiation; involved in regulation of synaptic plasticity [RGD, Feb 2006]