

## Product datasheet for **RR216153**

### Tnnt1 (NM\_001277262) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Tnnt1 (NM\_001277262) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Tnnt1  
**Synonyms:** Fang2; sTnT; tnTs  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR216153 representing NM\_001277262  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGTCAGACACCGAAGAACAAGAATATGAGGAGGAACAGGCAGAAGATGAGGAAGCAGTGGAGGAAGAAG  
AAGAAGAGCGCCCAAACCCAGCCGTCCTGTGGTTCCTCCTTTGATTCTCCGAAGATCCAGAAGGGGA  
GCGTGTGGACTTTGATGACATTCACCGAAGCGCATGGAGAAAGACTTACTGGAGCTGCAGACGCTCATC  
GATGTGCACTTTGAGCAGCGGAAAAAGGAAGAGGAGGAATCATTGCGCTAAAAGACCGTATTGAGAGGC  
GTCGCGCAGAGAGAGCTGAGCAACAGCGCTTCAGAACAGAAAAGGAGCGAGAACGCCAGGCCAAGCTGGC  
GGAGGAGAAGATGCGGAAGGAGGAGGAGGAGGCAAGAAGAGAGCAGAGGATGACGCCAAAAAGAAGAAG  
GTTCTGTCCAACATGGGGCTCATTTTGGGGGCTACCTGGTCAAGGCAGAACAGAAAGCGTGGGAAGCGAC  
AGACTGGTCGGGAGATGAAACTCAGGATCCTGTCTGAGAGGAAAAAGCCCTTGAACATCGACTACATGGG  
GGAGGACCAGCTCCGGGAGAAGGCCAGGAGCTGTCAGAATGGATCCACCAGCTGGAATCGGAGAAATTT  
GACCTGATGGAGAAATTAAGCAACAGAAATATGAGATCAATGTGCTGTACAACCGCATCAGTCATGCC  
AGAAATCCGAAAGGGGCTGGGAAGGGTCGAGTTGGAGGCCGCTGGAAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RR216153 representing NM\_001277262  
 Red=Cloning site Green=Tags(s)

MSDTEEQEYEEEAQAEDEEAVEEEEEERPKPSRPVVPPLIPPKIPEGERVDFDDIHRKRMEKDLLELQTLI  
 DVHFEQRKKEEELIALKDRIERRRAERAQQRFRTKERERQAKLAEEKMRKEEEEAKKRAEDDAK  
 VL SNMGAHFGGYLVKAEQKRGKRQTGREMCLRILSERKKPLNIDYMGEDQLREKAQELSEWIHQLESEKF  
 DLMEKLLKQKYEINLVYNRI SHAQKFRKGAGKGRVGRWK


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_001277262  
**ORF Size:** 750 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001277262.2</a> , <a href="#">NP_001264191.1</a>
<b>RefSeq Size:</b>	911 bp
<b>RefSeq ORF:</b>	753 bp
<b>Locus ID:</b>	171409
<b>UniProt ID:</b>	<a href="#">Q7TNB2</a>
<b>Cytogenetics:</b>	1q12
<b>MW:</b>	30 kDa
<b>Gene Summary:</b>	human homolog plays a role in the calcium sensitivity of contraction in skeletal muscle [RGD, Feb 2006]