

Product datasheet for **RR216154**

Tnnt3 (NM_001270672) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tnnt3 (NM_001270672) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Tnnt3
Synonyms: Tnt
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR216154 representing NM_001270672
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCTGACGAGGAACTGAACAAGTTGAGGAACAGTACGAAGAGGAAGAGGAGGCCAGGAGGAAGAGG
AAGTCCAAGAAGAGGAGAAAACCAAGACCCAACTTACTGCTCCTAAGATCCCGAAGGAGAGAAAAGTAGA
CTTCGATGATATCCAGAAGAAGCGTCAGAACAAGGACCTCATGGAGCTCCAGGCCCTATTGACAGCCAC
TTTGAAGCTAGAAAAGAAAAGAAGAGGAGCTAATTGCACTTAAGGAAAGAATTGAAAAGCGCCGTGCAG
AGAGGGCTGAGCAGCAAAGGATTCGCGCTGAGAAGGAGCGGGAGCGCCAGAACAGACTGGCGGAGGAGAA
GGCCAGAAGAGAGGAGGAAGATGCCAAGAGGAGAGCTGAAGATGACTTGAAGAAGAAAAGGCTCTGTCC
TCTATGGGTGCCAACTACAGCAGCTACCTGGCCAAGGCTGACCAGAAGAGAGGCAAGAAAACAGACAGCCC
GGGAAATGAAAAGAAGATTCTTGCTGAAAGGCGCAAACCTCTGAACATTGACCATCTTAGTGATGATAA
GCTGAGGGACAAGGCCAAGGAACCTCTGGGATACCTGTACCAACTGGAGACTGACAAATTTGAGTTGGG
GAGAAGCTGAAACGTGAGAAATACGATATTACCACCTCAGGAGCCGATTGACCAAGCCAGAAGCACA
GCAAGAAGGCCGGTGCCACAGCCAAGGGCAAAGTCGGCGGGCGCTGGAAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR216154 representing NM_001270672
 Red=Cloning site Green=Tags(s)

MSDEETEQVEEQYEEEEEAQEEEEVQEEKPRPKLTAPKIPEGEKVDFDDIQKKRQNKDLMELQALIDSH
 FEARKKEEELIALKERIEKRAERAQQRIRAEKERERQNLAEKARREEDAKRRAEDDLKKKALS
 SMGANYSYLAKADQKRGKKQTAREMKKKILAERRKPLNIDHLSDDKLRDKAKELWDTLYQLETDKFEFG
 EKLRQKYDITTLRSRIDQAQKHSKKAGATAKGVGGRWK


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001270672
ORF Size: 750 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001270672.1 , NP_001257601.1
RefSeq Size:	1053 bp
RefSeq ORF:	753 bp
Locus ID:	24838
Cytogenetics:	1q41
MW:	29.7 kDa
Gene Summary:	plays a role in muscle contraction [RGD, Feb 2006]