

Product datasheet for **RR216162**

Tnnt3 (NM_001270673) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tnnt3 (NM_001270673) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Tnnt3
Synonyms: Tnt
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR216162 representing NM_001270673
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCTGACGAGGAACTGAACAAGTTGAGGAACAGTACGAAGAGGAAGAGGAGGCCAGGAGGAAGAAG
TCCAGGAGGAAGCCCCGGAACCAGAGGAGAAACCAAGACCCAACTTACTGCTCCTAAGATCCCGGAAGG
AGAGAAAGTAGACTTCGATGATATCCAGAAGAAGCGTCAGAACAAGGACCTCATGGAGCTCCAGGCCCTC
ATTGACAGCCACTTTGAAGCTAGAAAGAAAGAAGAAGAGGAGCTAATTGCACTTAAGGAAAGAATTGAAA
AGCGCCGTGCAGAGAGGGCTGAGCAGCAAAGGATTCGCGCTGAGAAGGAGCGGGAGCGCCAGAACAGACT
GGCGGAGGAGAAGGCCAGAAGAGAGGAGGAAGATGCCAAGAGGAGAGCTGAAGATGACTTGAAGAAGAAA
AAGGCTCTGTCTCTATGGGTGCCAACTACAGCAGCTACCTGGCCAAGGCTGACCAGAAGAGAGGCAAGA
AACAGACAGCCCCGGAAATGAAAAGAAGATTCTTGTGAAAGGCGCAAACCTCTGAACATTGACCATCT
TAGTGATGATAAGCTGAGGGACAAGGCCAAGGAAGCTCTGGGATACCTTGTACCAACTGGAGACTGACAAA
TTTGAGTTTGGGAGAAGCTGAAACGTCAGAAAACGATATTACCACCCTCAGGAGCCGATTGACCAAG
CCCAGAAGCACAGCAAGAAGCCGGTGCCACAGCCAAGGCAAAGTCTGGCGGGCGCTGGAAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR216162 representing NM_001270673
 Red=Cloning site Green=Tags(s)

MSDEETEQVEEQYEEEEEAQEEVVQEEAPEPEEKPRPKLTAPKIPERGEKVFDDIQKKRQNKDLMELQAL
 IDSHFEARKKEEEELIALKERIEKRAERAQQRIAEKERERQNRLAEEKARREEDAKRRAEDDLKKK
 KALSSMGANYSSYLAKADQKRGGKQTAREMKKKILAERRKPLNIDHLSDDKLRDKAKELWDTLYQLETDK
 FEFGEKLRQKYDITTLRSRIDQAQKHSKKAGATAKGKVGGRWK


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001270673
ORF Size: 762 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001270673.1 , NP_001257602.1
RefSeq Size:	1065 bp
RefSeq ORF:	765 bp
Locus ID:	24838
UniProt ID:	P09739
Cytogenetics:	1q41
MW:	30.1 kDa
Gene Summary:	plays a role in muscle contraction [RGD, Feb 2006]