

Product datasheet for **RR216193**

Icam4 (NM_001172077) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Icam4 (NM_001172077) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Icam4
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR216193 representing NM_001172077
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGTCTGCCCTTCTGCTCCTTTGCTTTTGGTGGTGCCTATTTGAGAGGTGGAAGCTCCCAGC
ATGAGTGGATGCAAAGCCCTCCCGCACCTCCGTGACCTCAGTACCTTTCTTGGTGCCTTAGTCCGGA
GCTGGAGGCTGTGCCTCCCGGGGCTCAGTGTGGCTTAAGTGCAGCCACAAGTCCCGCTGCCGGTGCAT
TCCAGCCTGCGCACTCAACTGCGCAGGAAAGATAGTCAATGGATCCGGCTGGGTATCTTACCAGCTAC
TGGATGTGAGGGCCTGGAATCCAAGGTGCGCTGCGTCACTTGCAGGAGAAACCCGAGAGGCCAC
CGCCAGGATCACTGCTTTCAAACGGCCAGAAGTGTGATCTTGGAGCCTCCAGTCCAGTGGCCACAAG
TACACTCTGAGATGCTATGTGACACACGTGTTCCAGTGGGGTCTTGGTGGTGGCCTGAGGAAAGGTG
GCCGGTGGTTTATCATGAAAGTCTGCAGCGCTTACCGGCTCAGATTTGGCTAATGTCACGTTGACCTA
CGTGATGCGGGCCGGACTCAACGACCTTGCCAACCACTACCTGCCACGCGCCTCAATCTCGACGGG
CTAGTGGTGCAGCAGCTCAGCACCTGTTATGCTGACGGTCTCGCTTTAAGCCAGCCTCCATATCCT
TGGCCTCTATCTCCATCGCAACCCTGGTGGGGATCTTCTGGCTGTGGGGCTGTCTACGTGCGCAAGTA
CCTGGCCGTGCAGACTCAATTA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR216193 representing NM_001172077
 Red=Cloning site Green=Tags(s)

MESALLLLSLLLVAAYLRGGSSQHEWMQSPAPSVTSPVFLVRLSPELEAVPPGGSVLNC SHNCPLPVH
 SSLRTQLRQGI VNGSGWVSYQLLDVRAWNSKVRVCVTCAGETREATARITAFKRPRSVILEPPVLVGHK
 YTLRCYVTHVFPVGFLLVSLRKGRVVYHESLQRF TGSDLANVTLYVMRAGLNDLWQPLTCHARLNLDG
 LVVRSSSAPVMLTVLALSPASISLASISIATLVGIFLAVGAVYVRKYLA VQTQL


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001172077

ORF Size: 792 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001172077.1 , NP_001165548.1
RefSeq Size:	845 bp
RefSeq ORF:	795 bp
Locus ID:	298702
Cytogenetics:	8q13
MW:	28.7 kDa