

Product datasheet for **RR216289**

Hnrnpc (NM_001277601) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Hnrnpc (NM_001277601) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Hnrnpc
Synonyms:	hnRNP C; Hnrpc
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR216289 representing NM_001277601 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGCTAGCAATGTTACCAACAAGACAGATCCTCGATCCATGAATCCCCTGTATTCATTGGGAATCTCA
ACACTCTGGTGGTCAAGAAGTCTGATGTGGAGGCCATCTTTCAAAGTATGGCAAATTTGGGGTTGCTC
TGTGCATAAGGGCTTTGCCTTTGTCCAGTATGTTAATGAAAGAAATGCCCGAGCTGCTGTAGCTGGAGAG
GATGGCAGAATGATTGCTGGCCAGGTTTTAGATATTAACCTGGCTGCAGAGCCAAAAGTGAACCGAGGAA
AAGCGGGTGTGAAACGATCTGCAGCGGAGATGTACGGTTCCTCATTGACTTGGACTATGACTTTCAACG
CGATTATTATGACAGGATGTACAGTTACCCAGCACGTGTTCTCCTCCTCCTCCCATTGCTCGAGCTGTG
GTGCCCTTCCAAACGCCAGCGTGTGTCGGGGAACACCTCACGAAGGGGCAAAGTGGATTCATTCAAAGA
GTGGACAACGGGGATCTTCTTCCAAATCTGAAAGTTGAAAGGTGATGACCTTCAGGCCATTAAGAAAGGA
GCTGACTCAGATAAAACAAAAGTGGATTCTCTGCTGAAAGCCTGGAAGAAATGAAAAAGAACAAGC
AAACAAGCAGACTTGTCTTCTCGTCCCAGTAGAAATGAAGAATGAAAAGTCTGAAGAGGAGCAGAGCA
GTGCCTCTGTGAAGAAAGATGAGACTAATGTGAAATGGAGTCTGAGGCAGGTGCAGATGACTCTGCTGA
GGAGGGCGACCTGCTGGATGATGACGATAATGAAGATCGGGGGGATGACCAGCTGGAGTTGAAGGATGAT
GAAAAAGAGCCTGAGGAAGGAGAAGATGACAGAGATAGCCCAATGGGGAGGATGACTCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR216289 representing NM_001277601
Red=Cloning site Green=Tags(s)

MASNVTNKTDPRSMNSRVFIGNLNLTLVVKKSDVEAIFSKYGKIVGCSVHKGFVQYVNERNARAAVAGE
 DGRMIAGQVLDINLAAEPKVNVRGKAGVKRSAAEMYGSSFDLDYDFQRDYDRMYSYPARVPPPPPIARAV
 VPSKRQRVSGNTRRGKSGFNKSGQRGSSSKSGKLGDDLQAIKKELTQIKQKVDLSLLESLEKIEKEQS
 KQADLSFSSPVEMKNEKSEEEQSSASVKKDETNVKMESEAGADDSAEEDLLDDDDNEDRGDDQLELKD
 EKEPEEGEDDRDSANGEDDS


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001277601

ORF Size: 900 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001277601.1 , NP_001264530.1
RefSeq Size:	1763 bp
RefSeq ORF:	903 bp
Locus ID:	290046
Cytogenetics:	15p14
MW:	33.1 kDa