

Product datasheet for **RR216323**

Nek6 (NM_001277232) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Nek6 (NM_001277232) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Nek6
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR216323 representing NM_001277232
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCAGGACAGCCAGCCACATGCCTCATGGAGGGAGCCAAACCACCTCTGCCACGTGCTGGGGCCTG
CACACCTCCGGACCCACAGAGGCACCCAAACACCCTGTCTTTTCGCTGCTCACTGGCAGACTTCCAAAT
CGAGAAGAAGATTGGCCGAGGACAGTTCAGTGAGGTCTACAAGCCACCTGCCTGTGGACAGGAAGACA
GTGGCTCTGAAGAAGGTGCAGATTTTGTAGATGATGGATGCCAAAGCCAGACAGGACTGTGTCAAAGAAA
TTGGTCTCCTAAAGCAACTGAACCATCCAACATCATTAAGTACCTGGACTCCTTCATCGAAGACAATGA
GCTGAACATTGTGCTGGAGTTGGCTGACGCAGGTGACCTCTCACAGATGATCAAGTATTTCAAGAAGCAG
AAGCGGCTTATCCCTGAGAGGACTGTGTGGAAGTACTTCGTGCAGCTGTGCAGCGCTGTGGAGCACATGC
ATTCCCGACGTGTGATGCACCGAGACATCAAGCCTGCCAACGTGTTTCATCACAGCTACGGGCATCGTGAA
GCTTGGTGACCTCGGCCGGCCGCTTCTTCAGCTCAGAGACCACTGCAGCCACTCGCTAGTGGGAACA
CCATACTATATGTCACCCGAGAGGATCCACGAGAATGGCTACAACCTCAAGTCAGACATCTGGTCCCTAG
GCTGCCTGCTGTATGAGATGGCAGCTCTTCAGAGCCCCTTCTATGGAGATAAGATGAATCTTTCTCCCT
CTGCCAGAAGATTGAGCAGTGTGACTACCCGCCTCTCCGGGGGAGCACTACTCTGAGAAGCTTCGGGAA
CTGGTCAGCATGTGCATCTACCCGACCCTAACCCAGCCTGACATTGAGTATGTGCACCAAGGTGGCTA
AACAGATGCACGTCTGGACATCCAGCACC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

ORF Size:	939 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001277232.1 , NP_001264161.1
RefSeq Size:	2264 bp
RefSeq ORF:	942 bp
Locus ID:	360161
UniProt ID:	P59895
Cytogenetics:	3q12
MW:	35.8 kDa
Gene Summary:	may act as a serine/threonine protein kinase [RGD, Feb 2006]