

Product datasheet for **RR216444**

Scamp3 (NM_031724) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Scamp3 (NM_031724) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Scamp3
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR216444 representing NM_031724
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTCAGAACAGAGACGGCGGGAATCCATTCCCTGACTCCGGCGAGCTTGACAATCCCTTTAGGACC
CAGCTGTGATCCAGCACCGACCCAGCCAGCAGTATGCCACGCTTGACGTCTACAACCCCTTTGAGAACC
AGAGCAGCCCCACCAGCCTATGAGCCTCCTGCCAGCCAGCCCTCTGCCGCCACCCTCAGTCCT
TCTGTACAGTCCCAAGAAAGCTCAGCCCGCAGAGCCCAAGAACTACGGCTCCTACAGCACTCAGGCTA
CAGCTGCAGCAGCCACCGCTGAGCTGCTGAAGAAGCAGGAGGAGCTCAACCGAAGGCAGAGGAGTTGGA
CCGCAGGGAGCGAGAGCTGCAGCATGTTGCCCTGGGTGGTGCAGGCACTCGACAGAACAAGTGGCCTCCC
CTGCCGCTTTTTGCCCTGTGCAGCCTTGTCTTTTCCAAGACATCTCTATGGAGATCCCCAAGAATTT
AGAAGACAGTGTCTACCATGTAACCTCTGGATGTGCAGCACTCTGGCTCTTCTCCTGAACTTCTTCGC
CTGCCCTGGCCAGGTTCTGTGTGGACACCGGCAGTGGGTGAGGCTTTGGGCTGTCTATGCTCTGGCTGCT
CTTTTTACTCCCTGCTCCTTTGTCTGCTGGTACCGCCGATGTATAAGGCTTTCCGGAGTGATGTTTCAT
TCAATTTCTTCGTTTTTTCTTCATTTTCTTCGTCAGGATGTGTCTTTGTCTCCAGGCCATTGGTAT
CCCAGGTTGGGGCTTCAGTGGCTGGGTGTCCGCTCTCGTGGTGGTGGGTCACCCAGCCGTGGCGGTG
CTCATGCTGCTCGCCTGCTCTTCACTGGCATTGCCGTGCTGGGAATTGTGATGCTGAAGCGGATCC
ACTCCTTATACCGCCAAACAGGTGCCAGCTTTCAGAAGGCCAGCAAGAATTTGCTGCTGGTGTCTTCTC
CAACCTCGAGTTCGAAGTCCGCCCGCAATGCAGCTGCAGGGGCTGCTGAAAATGCCTTCAGGGCCCCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >RR216444 representing NM_031724
 Red=Cloning site Green=Tags(s)

MAQNRDGGNPFPPDSGELDNPFQDPAVIQHRPSQQYATLDVYNPFENREQPPPAYEPPAPAPLPPPSAP
 SVQSSRKLSPAEPKNYGSYSTQATAAAAATAELLKKQEELNRKAELDRRERELQHVALLGGAGTRQNNWPP
 LPSFCPVQPCFFQDISMEIPQEFQKTVSTMYLWMCSTLALLLNFFACLARFCVDTGSGSGFGLSMLWLL
 LFTPCSFVWCYRPMYKAFRSDSSFNFVFFVFFVFFVQDVFVFLQAIGIPGWGFSGWVSALVVVGSNPAVAV
 LMLLVALLFTGI AVLGI VMLKRIHSLYRQTGASFQKAQQEFAAGVFSNP AVR TAAANAAAGAAENAFRAP


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_031724

ORF Size:	1050 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_031724.1 , NP_113912.1
RefSeq Size:	1497 bp
RefSeq ORF:	1053 bp
Locus ID:	65169
Cytogenetics:	2q34
MW:	38.5 kDa
Gene Summary:	a putative membrane-trafficking protein that may be involved in endocytosis [RGD, Feb 2006]