

Product datasheet for **RR216452**

Otulin (NM_001302889) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Otulin (NM_001302889) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Otulin
Synonyms: Fam105b
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR216452 representing NM_001302889
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGAGTCGGGGAATCATGCCCAACCCGGAGCGTGCCCGGAGCGGGCTGCGCGGAGACGCCGGCGCGCG
AGGCAGGGGCTGCGGCGCGGGACGGCGGAAGGTGACGGTCGGCGCGCAGCCCCGGGCCGACCGGATG
CCCGGCCGAGCAGGAGGACATGTACCGTGTGCAGATGAAATAGAAAAGGAGAAAGAATTGCTAATA
CATGAAAGAGGGATATCAGCAGAACCCAGGTTAAGTGTGGCTCCTGAAATGGATATTATGGACTACTGCA
AAAAGGAATGGAGAGGAAATACCGAGAAAGCCATGTGCATGAGAAAGGGCTATGAGGAAGTTTCTCAGAA
GTTCACTTCCATAAGGCGAGTCCGTGGCGATAATTACTGTGCACTGCGGGCCACACTGTTCCAGGCTATG
AGCCAGCTGGCTGAGCTGCCCCCTGGCTGCAGGACCTGGAGCTCATACTGTTACCGAAAAGCTGATAA
ACAAATACCTGGATCAAGCAGTGAAACTTGGACTGAAATTTGATGGGAAGAGTGAGGACCTGGTTGA
AGAATCAAGGAATCCCTTGCCCTGCTAAGGAAGAAGTGGGTGAGCCTGGCTGCAATGGAAGCCGCTGAA
GCAAGGCAGACAGCCTGTGACGAGCTGTTCAAAACGAGGAGGAGGAGTACAGCCTCTACGAAGCTATGA
AGTTCCTGATGTTAGCCAGAGCCATCCAAGTGTATGATGACAAGGAAAAGGAAAGGAAGTTCCTTTTTT
CTCTGTGCTCTTGTGTTGCCGAGACACATCAAATGACCCTGAACAGCTCCTGAGGAACCACTGAACAG
GTGGGACACACTGGCGGCTTGAACAGGTTGAAATGTTCTTCTTGCCTATGCTGTTCCGCCACACCATCC
AGGTGTACCGGCTGTCCAAGTATAACACTGAGGAATTCATCACAGTCTATCCCACTGATCCCCCAAGGA
CTGGCCAATGGTGACCCTCATAGCTGAGGACGATCGGCACTATAACATCCCTGTCAGAGTGTGTGAGGAG
ACAAGTCTG

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR216452 representing NM_001302889
 Red=Cloning site Green=Tags(s)


MSRGIMPQPGAWPGAGCAETPAREAGAAARDGGKVTVGAQPRAATRCPAEHEEDMYRAADEIEKEKELLI
 HERGISAEPRLSVAPEMDIMDYCKKEWRGNTKAMCMRKYEEVSQKFTSIRRVRGDNYCALRATLFQAM
 SQLAELPPWLQDLELILLPEKLINKYTWIKQWKLGKFDGKSEDLVERIKESLALLRKKWVSLAAMEAAE
 ARQTACDELFTNEEEYSLEYAMKFLMLARAIQLYDDKEKGKEVPFFSVLLFARDTSNDPEQLLRNHLNQ
 VGHTGGLEQVEMFLLAYAVRHTIQVYRLSKYNTTEEFITVYPTDPPKDWPMTLIAEDDRHYNIPVRVCEE
 TSL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001302889

ORF Size:	1059 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001302889.1 , NP_001289818.1
RefSeq Size:	3311 bp
RefSeq ORF:	1062 bp
Locus ID:	100362554
Cytogenetics:	2q22
MW:	40.4 kDa