

Product datasheet for **RR216474**

Wars2 (NM_001168641) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Wars2 (NM_001168641) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Wars2
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR216474 representing NM_001168641 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGCTGTTCTCCGTGCGGAAGGCGCGTGAGTCTGGCGCTTCATCCGGGCATTCATAAAGGACCCA
CAGCTACTGCGGCTCCCCAGAAGGAGAGCGGGAGCGAGTGTCTCTGGCATTACAGCCTACAGGAATCCC
CCACCTGGGAAATTACCTTGGAGCCATCGAGAGCTGGGTGAAGTTACAGGAGGAGTATGACACGGTGATA
TACAGCATTGTGGACCTCACTCCATCACTGTCCCCAAGACCCCGGCATCCTGCAGCAGAGCATCCTGG
ACATGACTGCCGTGCTCCTTGCCTGCGGCATAGACCCAGAGAGAAGTATCCTTTCCAGCAGTCTCAGGT
GTCTGAACACACTCAGTTAAGTTGGATCCTCACCTGCATGGTGAGACTGCCTCGACTACAGCATTACAC
CAGTGGAAGGCGAAGGCTGCGAGGCAGAAGCATGATGGGACAGTAGGCCTGCTCACATACCTGTACTCC
AGGCGGCAGACATTCTGTGCTACAAGTCCACACACGTTCTGTGCGGGAAGATCAAGTCCAGCACATGGA
ACTAGTTCAGGATCTAGCTCGGAGTTTCAACCAGAAGTATGGCGAGCTCTTCCACTGCCAGGTCTATT
CTCACATCCATGAAGAAAGTAAATCTCTCGTGACCCTTCTGCCAAGATGTCAAAATCGGACCCTGACA
AACTGGCCACTGTCCAAATAACAGACAGCCAGAGGAGATTGTACGGAAATTCGCAAGGCTGTGACGGA
CTTCACGTCAGAAGTCACTACGAGCCGACAGCAGACCTGGTGTCTCCAACATGGTGGCGATCCATGCG
GCCGTGTAGGCCTCTCCGTGGAGGAGTGGTGCACAATAGTGCAGGCGTGGACTGCACGCTACAAGC
TGCTAGTGGCCGACGCTGTGATTGAGAAATTTGCTCCAATCAGAAGTAAATGAAAACTGAAATGGA
TAAGGACCATTTAAGAAAGGTTTACTTGTGGGATCTGCAAAAGCCAAAGAATTAGCCTCTCCTGTGTTCC
GAGGAGGTGAAGAAGTTGGTGGGATTCTG

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR216474 representing NM_001168641
 Red=Cloning site Green=Tags(s)

MALFSVRKARECWRFFIRALHKGPTATAAPQKESGERVFSGIQPTGIPHLGNYLGAIESWVKLQEEYDVI
 YSIVDLHSITVPQDPGILQQSILDMTAVLLACGIDPERSILFQQSQVSEHTQLSWILTCMVRLPRLQHLH
 QWKAKAARQKHGDTVGLLTYPVLQAADILCYKSTHVPVGEDQVQHMLVQDLARSFNQKYGELFPLPRSI
 LTSMKKVKSLRDPSAKMSKSDPKLATVQITDSPEEIVRKFRKAVTDFTEVITYEPDSRPGVSNMVAIHA
 AVSGLSVEEVVRNSAGVDTARYKLLVADAVIEKFAPIRSEIEKLMKDKHLRKVLLVGSAAKELASPVF
 EEVKKLVGIL


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001168641

ORF Size:	1080 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001168641.1 , NP_001162112.1
RefSeq Size:	1978 bp
RefSeq ORF:	1083 bp
Locus ID:	690654
Cytogenetics:	2q34
MW:	40.1 kDa