

Product datasheet for **RR216494**

Ctns (NM_001191647) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ctns (NM_001191647) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Ctns
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR216494 representing NM_001191647
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGATGAAGAAGTGGCTGCTTATTTTGACCTTTTCCTTCTGAAGTTCATAGAGAGAAGTGAATCAACCG
TCGGTCTCTGCTCCTCCACTGTGAAGCTGGTAAATGGCAGTTCAACCAACGTCAGCATTACCCTTGG
GCATCCACTAAATTCGACCTTGGTGATCACTTTTAAAATCACATTTAGTTCAAAAACCTTACTATTGTG
GAGCTTCCTGATGAGGTTATAGTGCCTCGTGGAGAGAAAAATGCCTCTTCCAAGTGACTTCTCAAATA
TTGGACAAGTGAAGTATCTGCATGGAATCATTCCAACAGACCTGCCCCAGAATCCAGTTCCTGGT
GATACACAGCAGAATTGTAGCGTCATAAACCAGGTGATTGGCTGGATCTACTTCATGGCCTGGTCCGTC
TCCTTCTACCCACAGGTGATCCAGAAGTGGAGACGGAAAAGTGTGCTGGTCTCAGCTTTGACTTCTTAG
CCCTGAAGTGGACAGGCTTCGTGGCCTACAGCGTTTTCAACATCGGTCTTTGTGGGTGCCCTACATCCA
GGAGGAGTTTCTCCTCAAATACCCCAATGGTGTGAACCTGTGACATTAACGACGCCTTCTTCAGCTTA
CATGCAGTCGCCCTCACCTGATTGTATCCTGCAGTGTGCCTGTATGAGCGAGGCAACCGCGGTGT
CACGGCCCTCCATTGGCTTCTGGTACTTGCCTGGCTTTCATCTTAGTCACCATGATTGTGGCTGCAGT
CGGTGTACCCACATGGCTGCAGTTCCTTCTGCTTCTCCTACATCAAGCTCATCATCACTGATCAAG
TATTTCCCCAGGCCATCATGAAGTACTTACTACAAAAGCACAAGGGTTGGAGCATTGGTGGTGTGCTCC
TGGATTTTACAGGGGGCAGCTTTAGCCTCCTCCAGATGTTCTCCAGTCTTACAATAATGACCAATGGAC
ATTGATCTTCGGAGATCCAACCAAGTTCGGACTTGGCGTCTTACCATCTTCTTTGACATTGTCTTCTTC
ATCCAGCACTTCTACTTGTACAGAAAGAAACCAGGGTACGACCAACTGGAC

ACGCGTACGCGGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR216494 representing NM_001191647
 Red=Cloning site Green=Tags(s)

MMKNWLLILTLFLLKFIERSESTVGLSAPPTVKLVNGSSTNVISITLGHPLNSTLVITFEITFSSKNLTIV
 ELPDEVI VPRGEKNASFQVTSQNI GQVTVYLHGNHSNQTCPRIQFLVIHSRIVSVINQVIGWIYFMAWSV
 SFYPQVIQNWRRKSVVGLSFDLALNLTGFVAVSVFNI GLLWVPYIQEEFLLKYPNGVNPVDINDAFFSL
 HAVALTLIVILQCCLYERGNQRVSRPSIGFLVLAWLFILVTMIVA AVGVTTWLQFLFCFSYIKLIITLIK
 YFPQAYMNFYKSTKGWSIGGVLLDFTGGSFLLQMFLQSYNNDQWTLIFGDPTKFGLVFTIFFDIVVF
 IQHFYL YRKKPGYDQLD


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001191647

ORF Size:	1101 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001191647.1 , NP_001178576.1
RefSeq Size:	2010 bp
RefSeq ORF:	1104 bp
Locus ID:	287478
Cytogenetics:	10q24
MW:	41.9 kDa