

Product datasheet for **RR216540**

Mbnl1 (NM_001191566) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Mbnl1 (NM_001191566) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Mbnl1
Synonyms: Mbnl
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR216540 representing NM_001191566
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGCTGTTAGTGTACACCAATTGGGACACAAAATGGCTAACACTGGAAGTATGTAGAGAGTTTCAA
GGGGACTTGCTCAGACCAGACACGGAATGTAATTTGCACATCCTTCGAAAAGCTGCCAAGTTGAAA
TGGACGAGTAATCGCCTGCTTTGATTCACTGAAAGGTCGTTGCTCCAGAGAAAAGTCAAATATCTTCAT
CCACCCCTCATTAAAAACACAATTGGAGATAAATGGGCGGAATAACTTGATTCAGCAGAAGAACATGG
CCATGCTGGCCCAGCAAATGCAGTTAGCCAACGCCATGATGCCCGGTGCCCACTGCAGCCCCTGCCAAT
GTTTTAGTTGCACCAAGCTTAGCCACCAATGCATCAGCAGCCTTTAATCCCTACCTGGGGCCTGTTTC
CCAAGCCTGGTTCCAGCAGAGATCTTGGCGACTGCACCAATGTTGGTCACGGGGAATCCTGGAGTTCCGG
TGCCAGCAGCTGCTGCAGCTGCTGCACAGAAGTTAATGCGGACAGACAGACTGGAGGTGTCCGAGAGTA
CCAGCGTGGCAATTGCAACCGAGGAGAAAATGACTGTGATTTGCTCATCCTGCTGACAGCACAATGATT
GATACCAATGACAACACAGTCACTGTCTGCATGGATTACATCAAGGGGAGATGCTCTCGGAAAAGTGCA
AATACTCCATCCCCCGCACACCTGCAAGCCAAGATCAAGGCTGCTCAATACCAGGTCAACCAGGCTGC
AGCAGCGCAGGCTGCAGCTACTGCAGCTGCCATGACTCAGTCGGCTGTCAAATCACTGAAGCGACCCCTC
GAGGCAACCTTTGACCTGGGAATTCCTCAAGCTGTACTTCCCCCATTGCCAAAAGAGGCTGCTCTTGAAA
AAACCAACGGTGCCACCGCAGTCTTTAACACGGGTATTTTCCAATACCAACAGGCTCTAGCCAACATGCA
GTTACAGCAGCATACAGCGTTCTCCACCAGTCCCATGGTGCACGGTGTACGCCAGCCACTGTGTCC
GCAGCAACAACATCTGCCACAAGTGTCCCTTCGCTGCAACAGCCACAGCCAACCAGATACCCATAATAT
CTGCCGAACATCTGACTAGCCACAAGTATGTTACCCAGATG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR216540 representing NM_001191566
 Red=Cloning site Green=Tags(s)


MAVSVTPIRDTKWL TLEVCREFQRTGCSRPDTECKFAHPSKSCQVENGRVIACFDSLKGRCsRENCKYLH
 PPPHLKTQLE INGRNLIQQKNMAMLAQQMQLANAMMPGAPLQPVPMFSVAPSLATNASAAFNPYLG PVS
 PSLVPAEILPTAPMLVTGNPGVPVPA AAAAAAQLMRTDRLEVCREYQRGNCNRGENDCRFAHPADSTMI
 DTNDNTVTVCMDYIKGRCSREKCKYFHPPAHLQAKIKAAQYQVNQAAAAQAAATAAMTQSAVKSLKRPL
 EATFDLGI PQAVLPPLPKRPALEKTN GATAVFNTGIFQYQQALANMQLQQTAF LPPVPMVHGATPATVS
 AATTSATSVPFAATATANQIPIISA EHLTSHKYVTQM

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001191566

ORF Size:	1161 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001191566.1 , NP_001178495.1
RefSeq Size:	1164 bp
RefSeq ORF:	1164 bp
Locus ID:	282635
Cytogenetics:	2q31
MW:	41.7 kDa