

Product datasheet for **RR216592**

Traf1 (NM_001271240) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Traf1 (NM_001271240) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Traf1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR216592 representing NM_001271240
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCCTCCAGCTCAGCCCTGATGAAAACGAGTTTCAATTTGGTTGTCCTCCCTGCTCCCAGCCAGGACC
CATCGGAGCCAGAGCTCTGCTGCACAGCCTGTCTCTGAGAACCTGAGAGACGATGAGGAACGGAT
CTGTCTAAATGCAGAGCAGACGGCCTCCAACCTGTGAGCCAGGAAGCCTTCTGACTCAGGAGAAGTT
CTCCCTGATGTGGCTGAGGCTGGAATCGTGTGCCCTTTCAGAGGTGTGGCTGTTCTTCAAGGGGAGCC
CACAGTCCGTGCAGGAGCATGAGGCCACCTCGCAGTCCTCACACCTGTACCTGCTTCTGGGGTCTAAA
GGAGTGGAAATCCTCACCAGGTCCCAACCTGGGCTCTGCACCCATGGCACTGGAGCAGAACCTGTGGAG
CTGCAGCTTCAGGAAGCTGTGGAAGTGACAGGGGACCTGGAGGTAGACTGCTACCGGCACCTTGTGTG
AGAGCCAGGAAGAAGCTGGCTTTCAGCACCTCCTGAAGGAGAAGCTGCTGGCTCAGCTGGAGGAGAAGCT
GCGTGTGTTCCGGAACATAGTTGCTGTCTCAACAAGGAAGTGGAGGCTTCCACCTGGCGCTGGCCGCC
TCCATCCACCAGAGTCAGTTGGACCGAGAGCACGTCCTGAACTTGGAGCAGAGGGTGGTGGAGTTACAGC
AGACCCTGGCTCAAAAAGACCAGGTTCTGGGCAAACCTTGAACAGTCTGCGACTCATGGAGGAGGCGTC
CTTTGATGGTACTTTCTGTGGAAGATACCAATGTACCAAGCGGTGCCACGAATCAGTGTGTGGCAGG
ACTGTCAGCCTTCTCTCCAGCTTTCTACACTGCCAAGTACGGTTACAAGTTGTGCCTGCGCCTGTACC
TGAATGGGGATGGCTCAGGCAAGAAGACCCACCTGTCCCTGTTTATTGTCATCATGAGGGGAGAATATGA
CGCTCTGCTGCCCTGGCCTTCCGGAACAAGGTCACCTTTATGCTACTTGACCAGAAACACCGTGAGCAC
GCCATTGATGCCTTCCGGCCTGACTGAGCTCAGCCTCCTCCAGCGCCGAGAGTGAGACCAACGTGG
CCAGCGGCTGCCCGCTCTTCTCCCCCTCAGCAAGCTGCAGTCGCCCAAGCACGCTACGTAAGGATGA
TACAATGTTCTCAAATGCATTGTGGACACCAAGTGT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR216592 representing NM_001271240
 Red=Cloning site Green=Tags(s)

MASSAPDENEFQFGCPPAPSQDPSEPRALCCTACLSENLRDDEERICPKCRADGLQPVSPGSLLTQEKV
 LPDVAEAGIVCPFAGVGCSPQSVQEHEATSQSSHL YLLGLVKWKSSPGPNLGSAPMALEQNLSE
 LQLQEAVEVTGDLEVDCYRAPCCSEQEELALQHLLKEKLLAQLLEKLRVFANIVAVLNKEVEASHLALAA
 SIHQSQLDREHVLNLEQRVVELQQTLAQKDQVLGKLEHSLRMEASFDTGTFLWKITNVTKRCHEVCGR
 TVSLFSPAFTAKYGYKLCLRLYLNGDGSKKTHLSLFIIVMRGEYDALLPWPFRNKVTFMLLDQNNREH
 AIDAFRPDLSSASFQRPQSETNVASGCPLFFPLSKLQSPKHAYVKDDTMFLKICIVDTSA


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001271240

ORF Size:	1227 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001271240.1 , NP_001258169.1
RefSeq Size:	1911 bp
RefSeq ORF:	1230 bp
Locus ID:	687813
Cytogenetics:	3p11
MW:	45.3 kDa