

## Product datasheet for **RR216598**

### Endou (NM\_001191069) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Endou (NM\_001191069) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Endou  
**Synonyms:** Pp11r  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR216598 representing NM\_001191069  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGAAAGCCTGTGTCCCTTAAATAGTGGCTATGCTCTGTGGCCTAACCTGGGCTGAAACCTGGAGTCGT  
GCGCATCTCGCTGTAATGAGAGATTTAACCGGGATGCTGTCTGCCAGTGTGATCGCCGGTCCCCAGCA  
TGATGACTGCTGTGACGACTATGAACATCTGTGCACTGATGAAGAAGGCCCCCCAGAGCCAAAGGCACTC  
CTGGATCCGGAGGATGAGATCCCAGCCAGCTACCTGTACTCCTCACCAAGCTCCTGCCAGGGCCGCTGCC  
GGGAAACCTACGACAAACACCATCCTTGCCACTGCAATGACCGCTGCCTAGAGTTTGGGAACTGCTGTGA  
GGATTTTGGCAACCTGTGTGGTGACCATGAGGGTTTCTCCACAGCAGCGATGCCATCACCAAGGAGGAG  
CTTAAGAGTATCTCCGAGACCATCTACAGGGTAGACATCAATAAAGCCCAGAAGAAAGACATTGTCCTCA  
ATAGTCAAACCAGATCTCCCGGCAGAGACAGGAAACCAAGTGGATCGATGCCAGAGCCGCTCTTCAA  
GTATGTCAACGAGAAGCTCTTCTCAAGCCTACCTATGCGGCCTTATCAACCTACTCAACAATTACCAG  
CGGGCCACAGGCCACGGGGAGCACTTCAGTCCCCAACAGCTAGAAGAGCAGGTTGTTTTCTCAGGGAGG  
TCATGAAGACAGCTGTCATGAAAGAGCTTTACAGCTTTTCCGTCACCAAAACCGTTACGGCTCCCAACA  
GGAGTTCTGTGGAAGATTTGAAGAACATGTGGTTTGGGCTCTACTCAAGAGGCAATGAAGAAGGGGATTCTG  
AGTGGCTTTGAACATGTCTTCTCTGCGCAGGTGAAGAAGGGCAAGTCACTGGCTTCCATAACTGGATCC  
GCTTCTACCTGCAAGAGAAGGAGGGTTTCTTACTATTACAGTCACAACCTACGACGGGCTTGGGATTC  
CTACCCAGATGTCTGGCCATGCAGTTCAACTGGGATGGCTATTACAAGGAAGTGGGCTCAGTCTTCATT  
GGCAGCAGCCCTGAGTTCGAGTTTGGCCTTTACTCTCTGTGCTTCATCACCAGGCCAGGCAAAAAGTGCC  
ATTTGAGCCTGGGGGCTATCCTTTGGCTATCCAGACCTACACCTGGGACAAGACCCTACGGGAATGG  
CAAGAAGTACATTGCCACAGCCTATGTGGTGTCTTCTACCCAA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RR216598 representing NM\_001191069  
Red=Cloning site Green=Tags(s)

MKACVPLIVAMLCGLTWAGNLESCASRCNERFNRDVAVCQCDRRCQPQHDDCCDDYEHLCTDEEGPPEPKAL  
 LDPEDEIPASYLYSSPSSCQGRCRETYDKHHPCHCNDRCLEFGNCCEDFGNLCDGHEGFSSHSSDAITKEE  
 LKSISETIYRVDINKAQKKDIVLNSQNQISPAETGNQVDRCPPEPLFKYVNEKLFASKPTAAFINLLNNYQ  
 RATGHGEHFSAQQLLEEQVVFLREVMKTAVMKELYSFLRHQNYRGSQQEFVEDLKNMWFGLYSRGNEEGDS  
 SGFEHVFSGEVKKGKVTGFHNWIRFYLQEKELLDYSSHNYDGPWDSYPDVLAMQFNWDGYYKEVGSVFI  
 GSSPEFEFALYSLCFITRPGKKCHLSLGGYPLAIQTYTWDKTTYGNGKKIATAYVVSSTQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_001191069

<b>ORF Size:</b>	1233 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001191069.1</a> , <a href="#">NP_001177998.1</a>
<b>RefSeq Size:</b>	1236 bp
<b>RefSeq ORF:</b>	1236 bp
<b>Locus ID:</b>	680317
<b>Cytogenetics:</b>	7q36
<b>MW:</b>	47 kDa