

Product datasheet for **RR216637**

Amdhd1 (NM_001191781) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Amdhd1 (NM_001191781) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Amdhd1
Synonyms: RGD1308879
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR216637 representing NM_001191781
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGC**C

ATGAATGGCCACCGCCTCCTGCTGGAGAACGCTCAGCAAGTGGTGTGGTGTGTGCCCGGGAGAGCTAT
 ATCTGACAGGGTCTGAGCTGCGCAGCCTGGCAGTGTGGAGGGCGCCAGCCTGGTGGTGGCCACTGATGG
 GCTTATAAAAGCTGTTGGCCCTGCTGACATTATCCGAAGACAGTTTTCCGGAGAAACGTTTGAAGAAAGA
 ATTGATTGCTCTGGAAAATGCATCTTGCCAGGTTTGGTGGACGGCACACACATCCGGTGTGGGCTGGTG
 AAAGAGTCCATGAGTTTGAATGAAGTTGGCAGGGGCCACGTACATGGACATCCACCAGGCTGGAGGAGG
 GATCAACTTACGGTGGAGCACACGCGCAAGCCTCGGAGGAGGAACTGTTCCGCTCCTTCCAGCAGCGG
 CTGCACTGATGATGAGGGCCGGCACCACACTGGTGGAGTGAAGAGCGGATATGGCCTCAACCTGGAGA
 CAGAGCTCAAGATGCTGCGTGTAAATCGAGCGTGGCCGGCGGGAGCTGCACCTCAGCCTCTCTGCCACTTA
 CTGTGGGGCTCACTCAGTACCCAAAGGGAAAACGCTGTGGAAGCTGCTGATGACATCATCAGTACCAC
 CTGCCAAGGCTGAAGGAACTCAGCGAAAAGGGGAACTGCATGTGGACAACATAGACGTTTCTGTGAGA
 AGGGCGTCTTTGATCTGGATACCACAGAAGGATTCTGGAAGGTGGAAAAAGGTGGGTTACAGATCAA
 CTTCCACGGGGATGAACTCCACCAATGAAGGCCGCTGAGCTTGGGGCGGAACTAGGAGCACGGGCGATC
 AGCCACCTGGAGGAAGTCAGCGATGAAGGCATCGCTGCCATGGCTGCGGCCAGTGTTCGCTGTCTTTC
 TGCCACCACGGCCTACATGCTGAGGCTGAAGCAGCCTCGCGCCAGGAAGATGCTAGATGAGGGAGTCAT
 AGTTGCTCTGGGAGCGATTTCAACCCCAACCGTACTGCTTTTCCATGCCTATGGTCATGCACCTGGCC
 TGTGTGAACATGAGGATGTCCATGCCAGAGGCTTTGGCCGCTGCCACCATCAATGCGGCTTATGCACTGG
 GAAAGTCTCACACCCACGGATCTTTGGAAGTTGGCAAGCAGGGGGATGCCATCATCAATGCATCCAG
 ATGGGAGCATTGATTTACCAGTTTGGAGGCCATCAAGAATTAATTGATTATGTTATAACTAAAGGAAAA
 GTCATCTATAAAAAG

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA



[View online >](#)


Protein Sequence: >RR216637 representing NM_001191781
 Red=Cloning site Green=Tags(s)

MNGHRLLENAQQVVLVCARGELYLTGSELRSLAVLEGASLVVGTDLIKAVGPADIIRRFQSGETFEER
 IDCSGKCILPGLVDAHTHPVWAGERVHEFAMKLAGATYMDIHQAGGGINFVETHQASEEELFRSFQQR
 LQCMRAGTTLVECKSGYGLNLETELKMLRVIERARRELHLSLSATYCGAHSVPKGTAVEAADDIISHH
 LPRLKELSEKELHVDNIDVFCEKGVFDLDTTRRILEGGKKVGLQINFHGDELHPMKAELGAELGARAI
 SHLEEVSDEGIAAMAAAQCSAVLLPTTAYMLRLKQPRARKMLDEGVIVALGSDFNPNAYCF SMPVMHLA
 CVMNRMSPEALAAATINAAYALGKSHTHGSLEVKGQGDAAIINASRWEHLIYQFGGHQELIDYVITKGG
 VIYKK


SGP TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-RsrII

Cloning Scheme:



Plasmid Map:



ACCN: NM_001191781

ORF Size:	1275 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001191781.1 , NP_001178710.1
RefSeq Size:	1290 bp
RefSeq ORF:	1278 bp
Locus ID:	299735
Cytogenetics:	7q13
MW:	46.6 kDa