

Product datasheet for **RR216688**

Efemp2 (NM_001277341) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Efemp2 (NM_001277341) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Efemp2
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR216688 representing NM_001277341
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGATCGCC

ATGCTCCCATTGCTCCTGCCTCCCGGGTCTTTACTGCTCTGGGCGTTTCTGCTGTTGCTCTTGGGAG
 CAGCGTCCCACAGGATCCGAAGAGCCGACAGCTACACGGAATGCACAGATGGCTATGAGTGGGATGC
 AGACAGCCAGCATTGCCGGGATGTCAATGAGTGCCTGACCATCCCCGAGGCTTGAAGGGTGAAATGAAA
 TGCATCAACCACTACGGGGCTATTTGTGTCTGCCTCGCTCTGCTGCCGTCATCAGTGATCTCCACAGTG
 AAGGACCTCCACCACCAGCAGCCTCTGCTCAACACCCAAATCCTTGCCCTCCGGGCTATGAACCTGATGA
 ACAGGAGAGCTGTGTGGATGTGGATGAGTGTACCCAGGCTATACATGACTGTGCGCCTAGCCAGGACTGC
 CATAACCTTCCTGGCTCCTATCAGTGCACCTGCCCGATGGTTACCGAAAAATTGGGCCGAATGTGTAG
 ACATAGATGAGTGTGCTACCGCTACTGCCAGCACCGATGCGTGAACCTGCCAGGCTCTTCCGATGCCA
 GTGTGAGCCCGCTTCCAGTTGGGACCTAACACCCGCTCTTGTGTGGATGTGAACGAGTGTGACATGGGA
 GCCCGTGTGAGCAGCGCTGCTCAACTCCTATGGGACCTTCTGTGCTGCTGTAACCAGGGCTATGAGC
 TACACCGGATGGCTTCTCCTGCAGTGATATCGATGAGTGCAGCTACTCCAGTTACTCTGCCAGTACCG
 CTGTGTCAATGAGCCAGTCTGCTTCTCCTGCTCACTGCCCAAGGCTACCAACTGCTGGCCACAAGGCTC
 TGCCAAGACATTGACGAGTGCAGAACAGGTGCACACCAAGTGTTCGGAGGCCAAACCTGCGTCAACTTCC
 ATGGGGTTACCCTGTGTGGACACCAACCGTTGTGTGGAGCCCTATGTCCAAGTGTGGACAACCGCTG
 CCTCTGCTCCTGCCTCAACCCCTTTGTGCGAGAGCAACCGTCAATCCATCGTGCACCGCTACATGAGCATC
 ACCTCAGAGCGAAGTGTGCTGCCGACGTGTTTTCAGATCCAGGCAACCTCTGTCTACCCTGGCGCTACA
 ATGCCTTTCAGATCCGTGCTGGGAACACACAGGGGACTTCTACATTAGGCAATCAACAATGTCAGCGC
 CATGCTGGTCTTGCAGACCGGTGACAGGACCCGGGAGTACGCTCTGGACCTCGAGATGGTACCATG
 AACTCCCTATGAGTTACCGGCCAGCTCTGACTGAGACTCACCGTCTTTGTGGGAGCCTATACCTTC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR216688 representing NM_001277341
Red=Cloning site Green=Tags(s)


MLPFASCLPGSLLLWAFLLLLLGAASPQDSEEPDSYTECTDGYEWDADSQHCRDVNECLTIPEACKGEMK
 CINHYGGYLCLPRSAAVISDLHSEGGPPPAASAQHNPNCPPGYEPDEQE SCVDVDECTQAIHDCRPSQDC
 HNLPGSYQCTCPDGYRKIGPECVDIDECRYRQCQHRVCNLPGSFRCQCEPGFQLGPNNRSCVDVNECDMG
 APCEQRCFNSYGTFLCRCNQGYELHRDGFSCSDIDECSYSSYLQCYRCVNEPGRFSCHCPQGYQLLATRL
 CQDIDECETGAHQCEAQT CVNFHGGYRCVDTNRCVPEYVQVSDNRCLCPASNPLCREQPSIVHRYMSI
 TSERSVPADVFIQATSYPGAYNAFQIRAGNTQGDIFYIRQINNVSAMLVLARPVTGPREYVLDLEMVTM
 NSLMSYRASSVLRLTVFVGAYTF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:	NM_001277341
ORF Size:	1329 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001277341.1 , NP_001264270.1
RefSeq Size:	1573 bp
RefSeq ORF:	1332 bp
Locus ID:	293677
Cytogenetics:	1q43
MW:	49.4 kDa