

## Product datasheet for **RR216739**

### Adcyap1r1 (NM\_001270582) Rat Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Adcyap1r1 (NM_001270582) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Adcyap1r1
Synonyms:	PAC1-R; PACAP-R1; PACAP-R1A; PACAPR1; PACAPR1A
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

**ORF Nucleotide Sequence:**

>RR216739 representing NM\_001270582  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGCATCGCC**

ATGGCCAGAGTCTCTGAGCTCTCCCTGACTGCTCTCTGCTGCCTGTGGCTATTGCTATGCACTCTGACT  
 GCATCTTCAAGAAGGAGCAAGCCATGTGCCTGGAGAGGATCCAGAGGGCCAACGACCTGATGGGACTAAA  
 CGAGTCTTCCCCAGGTTGCCCTGGCATGTGGGACAATATCACATGTTGGAAGCCAGCTCAAGTAGGTGAG  
 ATGGTCTTGTAAAGCTGCCCTGAGGTCTCCGGATCTTCAACCCGGACCAAGTCTGGATGACAGAAACCA  
 TAGGAGATTCTGGTTTTGCCGATAGTAATTCCTGGAGATCACAGACATGGGGTCTGGGCCGGAAGTCT  
 CACAGAGGACGGCTGGTGGAGCCCTTCCCCACTACTTCGATGCTTGTGGGTTTGTGATTATGAGCCT  
 GAGTCTGGAGATCAGGATTACTACTGTCCGTGAAGGCTCTACACAGTCCGCTACAGCACTTCCC  
 TCGCCACCCTCACTACTGCCATGGTTCATCTTGTCCGCTTCCGGAAGCTGCATTGCACTCGCAACTTCAT  
 CCACATGAACCTGTTTGTATCCTTCATGCTGAGGGCTATCTCCGCTTTCATCAAGGACTGGATCTGTAC  
 GCCGAGCAGGACAGCAGTCACTGCTTCCGTTCCACCGTGGAGTGCAAAGCTGTATGGTTTTCTTCCACT  
 ACTGCGTGGTGTCCAACACTTCTGGCTGTTTCATTGAAGGCTGTACCTCTTACACTGCTGGTGGAGAC  
 CTTCTTCCCTGAGAGGAGATATTTCTACTGGTACACCATCATCGGCTGGGGGACACCTACTGTGTGTGA  
 ACAGTGTGGGCTGTGCTGAGGCTCTATTTTGTATGATGCAGGATGCTGGGATATGAATGACAGCACAGCTC  
 TGTGGTGGGTGATCAAAGGCCCGTGGTTGGCTCTATAATGGTAACTTTGTGCTTTTCATCGGCATCAT  
 CATCATCTTGTACAGAAGCTGCAGTCCCAGACATGGGAGGCAACGAGTCCAGCATCTACTTACGGCTG  
 GCCCGCTCACCCCTACTGCTCATCCACTCTTCGGAATCCACTACACAGTATTGCGCTTCTCTCCAGAGA  
 ACGTCAGCAAGAGGAAAGACTTGTGTTGAGCTTGGGCTGGGCTCCTTCCAGGGCTTTGTGGTGGCTGT  
 ACTCTACTGCTTCTTGAATGGGAGGTACAGGCAGAGATTAAGAGGAAATGGAGGAGCTGGAAAGGTGAAC  
 CGTTACTTCACTATGGACTTCAAGCACCCGACCCGCTCCCTGGCCAGCAGTGGAGTAAATGGGGAAACC  
 AGCTGTCCATCCTGAGCAAGAGCAGCTCCAGCTCCGATGTCCAGCCTCCCGCCGACAACCTGGCCAC  
 C

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:**

>RR216739 representing NM\_001270582  
 Red=Cloning site Green=Tags(s)


MARVLQLSLTALLLPVAIAMHSDCIFKKEQAMCLERIQRANDLMGLNESSPGCPGMWDNITCWKPAQVGE  
 MVLVSCPEVFRIFNPDQVWMTETIGDSGFADSNLEITDMGVVGRNCTEDGWSEPFPHYFDACGFDDYEP  
 ESGDQDYIYLVKALYTVGYSTSLATLTTAMVILCRFRKLHCTRNFIMHNLVVSFMLRAISVFIKWILY  
 AEQDSSHCFVSTVECKAVMVFFHYCVVSNYFWLFI EGLYLFLLVETFFPERRYFYWYTIIGWGTPTVCV  
 TVWAVLRLYFDDAGCWMNDSTALWWVIKGPVVGSIIMVNFVLFIGIIILVQKLQSPDMGGNESSIYLR  
 ARSTLLLIFLFIHYTVFAFSPENVSKRERLVFELGLGSFQGFVAVLYCFLNGEVQAEIKRKRWSWKVN  
 RYFTMDFKHRHPSLASSGVNNGGTQLSILSKSSQLRMSSLPADNLAT

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**

**Plasmid Map:**


**ACCN:** NM\_001270582

**ORF Size:** 1401 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001270582.1](#), [NP\\_001257511.1](#)

**RefSeq Size:** 6151 bp

**RefSeq ORF:** 1404 bp

**Locus ID:** 24167

**UniProt ID:** [P32215](#)

**Cytogenetics:** 4q24

**MW:** 53.2 kDa

**Gene Summary:** regulates neural precursor proliferation; mediates inhibitory signaling for Shh-induced cerebellar granule precursor cell proliferation [RGD, Feb 2006]