

Product datasheet for **RR216867**

Cpne1 (NM_001256465) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cpne1 (NM_001256465) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Cpne1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RR216867 representing NM_001256465
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCCATTGTGTGACCTTGGTTCAGCTGTCCGTTTCCTGTGACCACCTCATTGACAAGGACATCGGCT
 CCAAGTCTGACCACTCTGCGTCTTTTACAGGATGTGGGAGGTGCCTGGGCTGAGCTTTGTCGGACTGA
 ACGTGTTTCGCAACTGCTCGAGCCCTGCATTCTCTAAGACTCTGCAGATAGAATACCACTTTGAGACTGTC
 CAGAAGTTACGCTTTGGAATCTATGACATAGACAACAAGACACCTGAGCTGGGGGATGATGACTTTCTAG
 GGGGAGCTGAATGCTCCCTAGGCCAGATTGTGTCCAGCCAGACGCTGACTCTACCCTTGATGTTGAAGCC
 TGGGAACTGCTGGACGGGAACCACTTACAGTTTCAGCTCAGGAGTTAAAGGACAGCCGTGTAGTAACC
 ATGGAGGTGGAAGCTAGAAACCTAGATAAGAAGGACTTTCTAGGAAAATCTGATCCATTCTGGAGTTCT
 TCCGGCAAGGTGATGGGAAGTGGCACCTGGCATATAGAAGTGGGTAGTCAAGAACAACCTGAACCTAC
 GTGAAACGCTTCTCAGTCTCTCTCCAGATTTCTGTGGCGGGGACCTCAACACACCCATCCAGGTGCGC
 TGCTCAGACTATGACAGTGTGGTGCATGATCTCATCGGTACCTTTCACACTACCTTGGCCCACTGC
 AAGCAGTCCCAGGCTGAGTTTGAATGTATCCACCCTGAGAAGCAGCAGAGAAGAAAGCTACAAGAAGT
 TGGAACTGTCTGTGTCAAGACTTCCCGGTAGAGACAGAGTATTCCTTCTGGACTATGTGATGGGAGGC
 TGCCAAATTAACCTTACCCTGGGTGTGGACTTCACTGGCTCCAACGGTGACCCATCCTCACCTGACTCCC
 TGCACTATCTGAGCCCTACAGGAGTCAATGAGTATCTGACAGCCCTGTGGAGTGTGGGCAGCGTGGTTCA
 AGACTATGACTCAGACAAGCTGTTTCCAGCATTGGATTGGAGCCCAAGTACCCCTGATTGGCAGGTC
 TCACATGAGTTTGCCTTGAACCTCAACCCAGTAACCCATACTGTGCAGGCATCCAGGGTATCGTAGACG
 CCTACCGTCAAGCACTTCCCAAGTTCTGTCTCTATGGCCCCACCACTTGCACCCATCATCAACCACGT
 GGCCAGGTTTGCAGCTCAGGCTGCACAACAGAGGACTGCCTCGCAATACTTCGTACTGCTGCTCTTGACT
 GACGGCGTGTGACCGATGTGGAGGCCAGTGTAAAGGCTGTGGTGAAGCCTCAAAGCTGCCATGTGCGG
 TGATCATCGTGGGTGTGGGTGGTGTGACTTTGAGGTCATGGAACAGCTGGATGCTGATGGCGGCCACT
 TCGTACACGCTCTGGAGAAGCAGCTGCCGAGACATAGTGCAGTTTGTGCCTTATCGACGATTTCAAGT
 GCTCCACGGGAGACACTTGAATGACTGTGCTTGCAGAAGTGCCCACTCAGATGGTTTCTACTTCAAAG
 CCCAAGTTGGGCCCATTAAGACTCCAGCCCCAGCAAGGGCCCTGCACAGGCCCCACAGGTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RR216867 representing NM_001256465
 Red=Cloning site Green=Tags(s)


MAHCVTLVQLSVSCDHLIDKDIGSKSDPLCVLLQDVGGAWAELCRTERVRNCSSPAFSKTLQIEYHFETV
 QKLRFGIYDIDNKTPELGDGDFLGGAECSLQIVSSQTLTLPMLKPGKPAGRGITVSAQELKDSRVVT
 MEVEARNLDKKDFLGKSDPFLFFRQGDGKWHLAYRTEVVKNNLNPTWKRFVSVLQHFCCGDLNTPIQVR
 CSDYDSGSHDLIGTFHTTLAQLQAVPAEFECIHPEKQQRKKSYSKNSGTVCKTCRVETEYSFLDYVMGG
 CQINFTVGVDFGTGNGDPSSPSLHYLSPTGVNEYLTALWSVGSVVQDYSDKLFPAFGGAQVPPDWQV
 SHEFALNFNPNPYCAGIQGIVDAYRQALPQVRLYGPTNFAPIINHVARFAAQAQRTASQYFVLLLLT
 DGAVTDVEATCKAVVEASKLPMSVIIIVGVGGADFEVMEQLDADGGPLRTRSGEAAARDIVQFVPPYRRFQ
 APRETLAMTVLAEVPTQMVSFYFKAQGWAPLKTLPAPAKGPAQAPQV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



ACCN: NM_001256465

ORF Size: 1608 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001256465.1](#), [NP_001243394.1](#)

RefSeq Size: 1611 bp

RefSeq ORF: 1611 bp


Locus ID: 362249

UniProt ID: [D4A1R8](#)

Cytogenetics: 3q42

MW: 58.8 kDa

Product images:



Circular map for RR216867