

## Product datasheet for **RR217042**

### **Srrm3 (NM\_001277394) Rat Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Srrm3 (NM_001277394) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Srrm3
Synonyms:	RGD1307391
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide  
Sequence:

>RR217042 representing NM\_001277394  
 Red=Cloning site Blue=ORF Green=Tags(s)

CTATAGGGCGGCCGGGAATTCGTCTGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCCGGCGC  
 GCC

ATGTCCTCCACCGTGAACAACGGGGCTACCAGCATGCCAGCTCCTCCGACGCTGCCAACGGCTTCCCGC  
 AGCCCGGAGGCTCGTCCGGGTCCTGGCCGCGCGGGAGGAGAGCTGCGCGCCGCCGAGCCAGGCTTGGT  
 GAAGCGCGCGCACCCGGGAGATCCTGGACCACGAGCGCAAGCGGCGGTAGAGCTCAAGTGCATGGAGTTG  
 CAGGAGATGATGGAGGAGCAGGGTACTCAGAGGAGGAGATCCGGCAGAAGGTCGGGACCTTCGGGCAGA  
 TGCTGATGGAGAAGGAGGGAGTACTGACAAAGGAGGACCGGCCTGGGGCCACATCGTGGCAGAGACTCC  
 GCGACGAATGGAGGGCTTGGAGCCCGCCTGGAGTACCCACCCTTCGATGAGGACGATGGCCCGTTGAT  
 TGGCAGTCCCGCTCGCTGCTACCGGGGACCGGGGTACAGGACCAAGCACTGGTCCAGCAGCTCCG  
 CATCGCCCCCTCCTAAGAAGAAGAAGAAAAGAAAGGCAGTCATCGGAGGAGTCGAAAAAGAGGAGACT  
 AGAGTCAGAATGCAGCTGTGGGAGTGTCCCCCTGCGCAAGAAGAAGAGAAGTGTGAAGAAACACCGT  
 CGAGATAGGTCTGATTCTGGGTGCGGAGAAAGAGACGATACAGATCGAGAAGCCTCAAGAGCAAACGAA  
 AAGAAAAATCAAAGAGCGAAAGAGGCCTCATACTGAGTCTCCAGGCCGAAGTCCCACCACCATAGCAG  
 CGTCAGCTCCCACAGCCCTCAGTGTCTCACACTACAGCGATTCCGGATCACCCAAACAGGCTGAGCCCC  
 AAGCATCGAGATGACGGACGGAAGACAGGCAGCCAGCGATCCAGCGGGAGCCGGTTCGCTGTCCCCATCGG  
 GCGGCAGCGGATGGGGTTCGCCCCAGCAGAACGGAGGCAGCAGGCAGCGAAGTGGAGCGCACGGGGGCCG  
 CCCCAGCTCGGCCACAGCCCGCCGATAAGCCAGCTCGCCCCGGGCTGCGACAAGGCAGCCGCCGCG  
 CCCACGCCCGCCGCGTGGGAAGGACAGCCAGAGTCCCCGCTCCGCGCCGTCGTCTCAGGGCCCGGGG  
 GCCGCGCGGGGAGGGACGGCCAGGCGTCCGCGCCGGCGGAGGAGGAGGAGCGATCGCGTCTCCGGC  
 CAACGCGCCCCGCGCAGGGGGCCCGGCGTGC AAAGCCCGCGCCGCTCGGGGCTCGTCTCGCTCGCTC  
 AGTGGGGCACACTCCAGCAGCGACTCCGGTGGCGCGCCCCGGGACCCGGGCCGAGCCGGGCTCCGAGC  
 GAGGTCACGGCGGACACGAAAACGGGCAAGGAGCGACCCCGCGAGCACGGCCCGCAAGCACCTCGCC  
 ATCACCAGGCGCGCATGGCCGACGGGGCGGCCAGAGGGGAACAGTTCGTGCGCAGCCCGGGCCCCCAT  
 CCCGGGTCGTGGAGTCCAGCCGCTCGCCCTCAAATCCCGCTCGCGTCCCCAGACAAAAGGACCCGCA  
 GCCCAGCCTCTACCGTCTCCAAGAAACCTTTGGCCGGGACAAGACAGCGAAGGCCGTGCCAGGCA  
 TGCCGAGGCTGAGGCCGCCGACCCGGGCTCGTCCCGAGCTACTCGCCGATCCGAAAGCGGCGTCCG  
 GACTCGCCGAGTTTCATGGAGCCGAGACGCATCACAGTCCCGAAAGCGTCCCATCCCCTACTACCGGC  
 CCAGCCCCCTCTTCTCCAGCTGCTTGAGCAGCGACTACTCCAGCCGAGCCGAGCAGGAGCCCAAG  
 CCTGGGCACAGCCATGGAAGCTACAGCAGCCGAGTACGGGACTCGCAGCCGTTCTGACGCGCTCC  
 AGGAGCCGAGTCCCAGTACCACAGCCGGAGCAGCTCTGAGAGTGGGGGCTTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RR217042 representing NM\_001277394  
 Red=Cloning site Green=Tags(s)


MSSTVNNGATSMPPAPDAANGFPQPGGSSGSWPRAEEELRAAEPGLVKRAHREILDHERKRRVELKCMEL  
 QEMMEEQGYSEEEIRQKVGTFGQMLMEKEGVLTKEDRPGAHVAETPRRMEGLEPGLYPPFDEDDGPVD  
 CDCPVACYRHRGYRTKHWSSSASPPPKKKKKKGGSHRRSRKKRRLESECSCGASPLRKKKRSVKKHR  
 RDRSDSGSRRKRRYRSRLSKRKEKIKERKRPHTESPGRRSHHHSSVSSHSPSVSSHYSDSGSPNRLSP  
 KHRDDGRKTGSQRSSGSRLSPSGSGWGSPQQNGGSRQSRGSAHGGRPGSAHSPDPKSSPRACDKAAAA  
 PTPPARGKDSQSPRSAPSSQGRGGAAGGTARRRRRRRRRRSRSSANAPRRRRRRRAKPAPPRGSSRSL  
 SGAHSSSDSGGGAPGPGPEPGERGHGGHKKRAKERPPRARPASTSPSPGAHGRRGGPEGNSSSRSPGPH  
 PGWSSSRSPSKRSRSPDKRTRSPSLSPSPKKPLGRDKDSEGRARHAEAEAARTRRRSYSPIRKRRR  
 DSPSFMEPRRITSARKRPIPYRPSSSSSSCLSSDYSSRSRSPSPGHSHGSYSRSRSHGTRSRSCSAS  
 RSRSPSYHSRSESSEGGF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:**

AscI-MluI

**Cloning Scheme:**



**ACCN:** NM\_001277394

**ORF Size:** 1944 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001277394.1](#), [NP\\_001264323.1](#)

**RefSeq Size:** 3549 bp


**RefSeq ORF:** 1947 bp

**Locus ID:** 685890

**Cytogenetics:** 12q12

**MW:** 70.8 kDa

**Product images:**



Circular map for RR217042