

Product datasheet for **RR217050**

Ddx3 (NM_001167665) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ddx3 (NM_001167665) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ddx3
Synonyms:	Ddx3y
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RR217050 representing NM_001167665
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGCGCGGAAAATACGGTTGGCCTTGACCAGCAATTTGCTGGTCTTGACATGAACTCAGATAATCAGA
 ATGGAGGAGTTAGTACAGAAAGCAAAGGCGCTATATACTCCTCACTTAAGGAACAGAGAAACATCTAA
 AGGTGTCTGTGATAAAGGACAATTCAGGGTGGAGTAAAGATAAAGATGCCTACAGCAGTTTTGGATCTCGG
 GATTCAGAGGAAAACCAATTAATTCAGTGATCGTGGAGTGGATCCAGAGGAAGGTTTGATGATCCTG
 GTCGAAGTGACTATGATTGATTGCCAGTCGTGACCGAAGTGGATTTGGCAAATTTGAACGAAGTGGTTA
 TAGTCGTTGGGGTGACAAATCAGATGAAGATGATTGGTCAAACCACTCCACCAAGTGAACGCTTAGAG
 CAAGAATTTTCTCTGGAGGTAATGCAGGGATTAACCTTGAGAAATATGATGATATACCAGTAGAGGCAA
 CTGAAATAATGTCTCCTCACATTGAAAATTCAGTGATGTTGAGATGGGAGAAATTATCATGGGAAA
 CATTGAACTTACTCGTTACACTCGTCTACTCCAGTTCAGAAACATGCTATTCCAATTATTAAGAAAAA
 AGGGACTTAATGGCTTGTGCTCAAACAGGGTCTGGATGTACTGCTGCATTTCTACTACCCATTCTCAGTC
 AGATTTACACAGATGGTCCAGGAGAGGCTTTGAAGGCTATGAAGGAAAACGGAAGATATGGTCGCCGTAA
 ACAATACCCAATCTCCTTGGTGTAGCTCCAACAAGAGAATTGGCTGTGCAGATCTATGAGGAAGCCAGA
 AAATTTTCGTACCGGTCTAGAGTCCGACCTTGTGTAGTGTATGGTGGTGTGATACTGTTCAACAGATTC
 GGGACTTAGAACGAGGATGCCACTTGTAGTTGCCACACCAGGACGTCTAGTGGACATGATGAAAGAGG
 AAAGATTGGACTAGATTTCTGCAAACTTAGTGTGGATGAAGCTGACAGGATGCTGGATATGGGATTT
 GAACCTCAGATACGTCGTATAGTTGAACAGGACACTATGCCACCAAGGGAGTTCGTACACCATGATGT
 TTAGTGCTACTTTTCTAAGGAGATACAGATGCTTGTGCTGACTTTTTAGATGAATACATCTTTCTGGC
 TGTAGGCAGAGTAGGCTCTACTTCTGAGAACATCACACAGAAAGTAGTATGGGTTGAAGAGTTAGACAAA
 CGGTCATTTTTACTTGATCTCTTAAATGCAACAGGAGGGATTATTGACTTTAGTATTTGTGGAGACTA
 AAAAGGGTGCAGATTCCTGGAGAATCTTATTTCAAGAAAGATATGCATGACTAGTATTCATGGAGA
 CAGGTCACAGAAAGATCGAGAGGAGGCCCTTACCAATTCGCTCAGGGAGGAAGCCAATTTTAGTGGCT
 ACAGCTGTGGCAGCAAGAGGTCTAGACATTTCAAATGTGAAACATGTTATCAATTTTGTGGCAAGTG
 ATATTGAAGAATATGTGCATCGTATTGGCCGCACAGGACGCGTAGGAAACCTTGGTCTTGCCACCTCCTT
 TTTCAATGAGAGGAATTTGAACATTACGAAAGATTTGTTAGATCTCCTCGTTGAGGCTAAACAAGAAGTA
 CCATCTTGGTTGAAAGTATGGCTTATGAACACCACTACAAGGAAACAGCCGGGACGTTCTAAAAGCA
 GATTCAGTGGAGGATTTGGTGCCAGAGACTATCGACAGAGCAGTGGTTATGCAAATTCGGGTTAATAA
 TAATCGTGCCAACAGCAGCCGAAGTGGTAGTACCCACAACAGAGGATTTGGTGGAGGTGGCTATGGA
 GTTTTCTACAGCAATGAAGGATATGGAGGAAATTATAACTCCCAGGCAGTTGACTGGTGGGGCAAT

AG**GCGACCG**ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR217050 representing NM_001167665
Red=Cloning site Green=Tags(s)

```


MAAENTVGLDQQFAGLDMNSDNQNGGVSTESKGRYIPPHLRNRETSKGVCDKDNSGWSKDKDAYSSFGSR
DSRGKPNYFSDRGSGSRGRFDDPGRSDYDCIASRDRSGFGKFERSGYSRWGDKSDEDDWSKPLPPSERLE
QELFSGGNAGINFKEYDDIPVEATGNCCPPHIENFSDVEMGEIIMGNIELTRYTRPTPVQKHAIPIIKEK
RDLMACAQGTSGCTAAFLLPILSQIYTDGPGEALKAMKENGRYGRRKQYPI SLVLAPTRELAVQIYEEAR
KFSYRSRVRPCVYVGGADTVQQIRDLEGRCHLLVATPGRLVDMMERGKIGLDFCKYLVLEADRMLDMGF
EPQIRRIVEQDTMPPKGVRRHTMMFSATFPKEIQMLARDFLDEYIFLAVGRVGTSENITQKVVVVEELDK
RSFLDLLLNATGRDSLTLVVFVETKKGADSLNLFQERYACTSIHGDRSQKDREEALHQFRSGRKPI LVA
TAVAARGLDISNVKHVINFDLPDIEEYVHRIGRTGRVGNLGLAT SFFNERNLNI TKDLLLVEAKQEV
PSWLESMAYEHYKGNRGRSKSRFSGGF GARDYRQSSGYANS GFNRRANSRRTSGSTHNRGFGGGGYG
GFYSNEGYYGNYNSQAVDWWGN
    
```

SGPTRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-RsrII

Cloning Scheme:



ACCN: NM_001167665

ORF Size: 1956 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001167665.1](#), [NP_001161137.1](#)

RefSeq Size: 3087 bp


RefSeq ORF: 1959 bp

Locus ID: 100312982

Cytogenetics: Yq11

MW: 73 kDa

Product images:



Circular map for RR217050