

Product datasheet for **RR217066**

Pus7 (NM_001170589) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pus7 (NM_001170589) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Pus7
Synonyms:	RGD1307054
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RR217066 representing NM_001170589
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGAGATGACAAGTGTGTCCCTGAAACGCGGGTGCCTCGTTGTGGAAGATAACGACAGTGTGACCCAC
 ATGAGGAGACGAAGAAACAAGGTGTCTGAGGGCTGTGTGACCAACAGTCCGGACGGGGTGGAGAGTGA
 CAGTCTTCACAGCAGTGAGACCAAGCCGGGCCCCAGAAAGCGAAAGCACCTTGAAAGATGGCGAGAAC
 TCCGGTGTCTCAGGTGCAGGAGGAAGAGGAGGAGGAGGAAGAGGACGGCCTTTCGGAGGCAGGGGAAG
 AGGAGGAGCCAGAGAGCTTTCGGGACATGATGAAGCATGGGCTCACCGAAGTGTGTGGGCATCACCAA
 GTTCGTGAGTCTCATCAGGGCTTCTCGGGATCTTAAAAGAAAGATACTCTGACTTCGTGGTTCATGAA
 ATCGGGAAGGATGGGCGGATCAGCCATTTGGATGACTTATCTGTTCCAGTGGACGAAGAGGACCCTCCAG
 AAGATGCCCTCACTGTCTTGACAGCTGAGGACAGGCAGCAGTTGGAGGAGCTCCAGCTTTTAAAAATAA
 AGAAACAGTGTGCCATCGAGGTCAATGAGGACACCAAAAGAGAAACGAACCATCATCCATCAAGCCATT
 AAATCTCTCTCCAGGACTAGAGACAAAAACAGAGGAGAGAGGGGAAGAGATACATTGTAGCCTACC
 ATGCTGTGGCAAGAAGGCTTTGGCAAATCCAAGAAAACATTCTTGCCAAAATCCCGGGGAAGTACTG
 CCACTTTGTTCTCTACAAGGAAAACAAAGACACCATGGACGCCATCAATGTGTTGTCCAAGTATCTGCGG
 GTAAAACCAAACATATTCTCCTATATGGGAACCAAGACAAAAGGGCTATCACAGTCCAGGAGATCGCTG
 TTCTCAAATAAGTGCACAAAGGCTTGGCCACCTGAATAAGTGCCTGATGAACCTTAAGTTAGGCAATTT
 CAGCTATCAGAAAACCCCTCTGAAGTTAGGAGAGCTTCAAGGAAACCACTTCACTGTCTGTTCTCAGGAAT
 ATAACGGGAACGGATGAACAAGTGCAGCAAGCCATGCATTGCTGCAGCAGACCGGCTTCATCAACTACT
 ACGGCATGCAGAGGTTTTGGGACCACGGCTGTTCCCTACGTACCAGGTCGGAAGAGCAATTCACAGAATTC
 CTGGACAGAAGTCATGGATTTAATACTGAAACCCGCTCTGGAGCTGAGAAGGGTTACTTGGTGAATGC
 AGAGAAGAGTGGGCCAAGACCAAGACCCAGCCTCTGCCCTCAAAAAGCTACCCGTCAAAAGGTGTGTGG
 AAGGGCAGCTGCTTCGAGGACTTCAAAGTATGGACTGAAGAATATAGTCTCTGCGTTTGGCATAATACC
 AAGAAACAACCGTTAATGTACATCCACAGCTACCAGAGTTATGTGTGGAATAGCATGGTGAGCAGGAGA
 ATCGAGCTGTATGGGCTGAAGCCTGTGCCGGGGACCTCGTGCTCAAAGGAGCCACGGCCACCTATATTG
 AGGAAGATGACATTAATAACTACTCCATTCACGATGTGGTGTGCCTTTGCTGGCTTTGATGTCATCTA
 CCCCAAACATAAAATTAGCGAGGCTACAGGAAATGCTTGCCGCAGACAACCTTGATATTGACAACATG
 AGACACACGATTGAGATTACTCCTTATCGGGGCTACCGAAAGATCATCATCCGACCTCAGAGCGTCA
 GCTGGGAAGTGGTTGCATACGATGATCCAAAATTCACCTTTTCAACACGGATGTGGACAACCTCGAGGG
 GAAACCACCACAGTCTTTGCTTCTGAAGGCAGATACAGGGCTCTGAAGATGGACTTCTCTCTCTCC
 TCTACTTACGCCACCATGGCCATTGAGAAGTGCTCAAAATGGACACCAGCATCAAGAACCAGACCCAGC
 TAAACACGAGCTGGCTGCGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR217066 representing NM_001170589
Red=Cloning site Green=Tags(s)

```
MEMTSVSLKRGCLVVEDNDSVTPHEETKKQKVSEGCVTNSPDGVEDSLHSSETKPGPPEAESTLKDGEN
SGAQVQEEEEEEEEEDGLSEAGEEEPEFSADMMKHGL TELDVGITK FVSPHQGFSGILKERYSD FVVHE
IGKDGRISHLDDL SVPVDEEDPPEDAL TVLTAEDRQQL EELQLFKNKETSVAIEVIEDTKEKRTIIHQAI
KSLFPGLETKTEEREGKRYIVAYHAAGKKALANPRKHSWPKSRGSYCHFVLYKENKDTMDAINVLSKYLR
VKPNIFSYMGTDKRAITVQEIAVLKISAQRLAHLNKCLMNFKLGNF SYQKTPLKLGELQGNHFTVVLRN
ITGTDEQVQAMHSLQQTGF INYYGMQRF GTTAVPTYQVGRAILQNSWTEVMDLILKPRSGAEKGYLVKC
REEWAKTKDPASALKKLPVKRCVEGQLLRGLSKYGLKNI VSAFGIIPRNNRLMYIHSYQSYVWNSMVSRR
IELYGLKPVPGDLVLKGATATYIEEDDINNYSIHDVVMPLPGFDVIYPKHKISEAYREMLAADNLDIDNM
RHTIRDYSLSGAYRKIIIRPQSVSWEVVAYDDPKIPLFNTDVDNLE GKPPPVFASEGRYALKMDFSLPP
STYATMAIREVLKMDTSIKNQTLNTSWLR
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001170589

ORF Size: 1980 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001170589.1](#), [NP_001164060.1](#)

RefSeq Size: 2464 bp


RefSeq ORF: 1983 bp

Locus ID: 296751

Cytogenetics: 4q11

MW: 74.6 kDa

Product images:



Circular map for RR217066