

Product datasheet for **RR217069**

Mepce (NM_001100678) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Mepce (NM_001100678) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Mepce
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RR217069 representing NM_001100678
 Red=Cloning site Blue=ORF Green=Tags(s)

CTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCCGGCGC
 GCC

ATGATCGAGATGGCGCGGAGAAGGAGCCCTTTCTGGTGCCAGCCCCACCGCCCGCTCAAAGATGAGT
 CGGGTGGAGGGGCGGCCCGAGGTGCGATCGCACCAAGAGGCTGCCTCCGGGGAGCTCCTCGACGGGAC
 AGAACGCGGGCCGGGCCACGCGCGCACTCGGCCGGGGCCACAGCCAGCGGAGGGGGCGCCCCAGGCA
 CAGACGCATGGGAGCCCCACGGGCGAGCTGCCGCTCCTGCGGACGTGGGGGAGGAGCGCCGGGAGGGG
 GCGGGACGGACCTGGGTCCCCAGCCCCCTCTCGACCTCGCAATGGCTACCAGCCCCACCAGCCCCAGT
 GGGTGGTGGGGCAAAGAAGAAATAGCTGTAACGTGGGGGAGGGAGCGGTGGAAGCTTCAAACATCCG
 GCCTTCAAGAGGCGCCGGCGAGTGAATTCGACTGCGACTCCGTGTTACCATCCAATTTCTTCTGGGGG
 GTAATATATTTGATCCACTGAACCTGAATAGCCTCCTGGATGAGGAAGTGAGTCGCGCACTCAATGCGGA
 GACCCCAAAGTCATCCCCACTTCGGCCAAAGGGGAGAGACCCTGTGGAGATCCTCATCCCAAAGATATC
 ACTGACCCACTCAGTCTCAACACGTGCACTGATGAAGCCATGTAGTGCTTGCTTACCACCTCAAGATTG
 GTCGAAGCGCCATAGACACCGGGACCGCACCACCAGCAGCAGCAGGCATCTGGAGGGAATGATAGCCA
 AGCTGTGCTGCCACTGACCCTCTCACACCCTCAGTCCATGGGGAGGGTGTACGCAGCAGCACAGAGGC
 CAGAACCAGATGACCCTCAACCCTATGAACTCAACACAGCCATCAACTGCAGGGATGAAGTTGTGTCTC
 CCCTTCGCTGCCCCTGCAGGGTTCCCTCAGGCTCCCTTCTGCCCTCCAGCTGCCTCAGTTACCTCTGC
 ACCCTCTCTCTCTCTCCGACATCGCAAACGTGCGAGGACTCCAGCAAGTCAGAGGCAGGGGCTAGG
 GGTGGAAGCCAGGGTTCCAAGGAAAAGGGCAGAGGGAGTGGGGGAGGCCGTACCACCACCACCTACCTG
 CAACTGGCTTCAAAAAGCAACAGCTCAAGTCCAGTATGGGAATTACTGCAAGTACTATGGCTATCGCAA
 TCCTTCTGTGAGGACATGCGCCTTAGGGTATTGAAGCCTGAGTGGTTTCAAGGCCGGGACGTCTAGAT
 CTGGGCTGCAACGTTGGCCATCTGACCCTAAGTATTGCCTGCAAGTGGGGCCAGCCCGCATGGTGGGGC
 TGGACATTGATCCACGACTTATCCACTCTGCCCGCAAAATATCCGACACTACCTGTCCGAGGAGCTTCG
 TCTCAAGCTCAGACTTCTGAGGGGACCCAGGGACAGAGGGGAGGAAGGGACTACCACTGTCCGAAAA
 AGGAGCTGCTTCCAGCCTCACTGACAGCCAGCCGGGTCCCATCGCTGCACCCCAATTTGCCCTGGATG
 GAGCGGACACATCTGTCTCCCCAACAATGTCGTCTTCGTCACGGGTAATATGTCCTGGACAGAGATGA
 GCTGGTGGATGCCAAAGACCTGAGTATGATGTGGTCTCTGCTTCAGCCTCACCAAGTGGGTGCACCTG
 AACTGGGGTGTGAGGGGCTGAAGAGAATGTTCCGCCGATCTACCGGCATCTGCGACCTGGGGCATTCT
 TGGTCTTGAGGCCCAACCGTGGTCTGCTACTGCAGGAGAAAGTCTCTCACAGAAACAATCTATAAGAA
 CTACTTTCGAATCCAGCTGAAGCCAGAACAGTTACGCTCCTACCTGACATCTCCAGAGGTGGGTTTCTCC
 AGCTATGAGCTTGTGGCCACACCAAGTAACACCTCCAGAGGCTTCCAGCGTCCGGTGTACTTGTCCACA
 AGGCCGATCTCCAGCCAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR217069 representing NM_001100678
 Red=Cloning site Green=Tags(s)


MIEMAAEKEPFLVPAPPPLKDESGGGGPEVRSHQEAASGELLDGTERGPGPRAHSAGATASGGGGPQA
 QTHGEPHGAAAPADVGEERRGGGDTLGPAPRRPRNGYQPHRPPVGGGKRRNSCNSVGGSGGSFKHP
 AFKRRRVNSDCDSVLP SNFLLGGNIFDPLNLNSLLDEEVSRLNAETPKSSPLPAKGRDPVEILIPKDI
 TDPLSLNTCTDEAHVVLASPLKIGRKRHRHRGPHHQQQASGGNDSQAVLPTDPLTPSVHGEGATQQHRG
 QNRDDPQPYELNTAINCRDEVVSPALQSSGSL SAPPAA SVTSAPSSSSSRHRKRRRTSSKSEAGAR
 GGSQGSKEKGRSGGGRHHHPLPATGFKKQQLKFQYGNKYGYRNPSCEDMRLRVLKPWFQGRDVL D
 LGCNVGHLTLSIACKWGPARMVGLDIDPRLIHSARQNI RHYLSEELRLQAQTSEGDPGTEGEEGTTVVRK
 RSCFPASLTASRGPAAAPHLPLDGADTSVFPNNVVFVTGNVYLDRELVDAQRPEYDVVLCFSLTKVHHL
 NWGDEGLKRMFRRIYRHLRPGGILVLEPQPWSSYCRKSLTETIYKNYFRIQLKPEQFSSYLTSPEVGF S
 SYELVATPSNTRS RGFQRPVYLFHKARSPSH

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites:

AscI-MluI

Cloning Scheme:



ACCN: NM_001100678

ORF Size: 1980 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001100678.1](#), [NP_001094148.1](#)

RefSeq Size: 3120 bp


RefSeq ORF: 1983 bp

Locus ID: 304361

Cytogenetics: 12q12

MW: 71.6 kDa

Product images:



Circular map for RR217069