

Product datasheet for **RR217077**

Ttc30a (NM_001271494) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ttc30a (NM_001271494) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ttc30a
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RR217077 representing NM_001271494
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCGTGGCAAAGCAGCTCGAAAGTCCCGACGGTGAGTTCACAGCGGTGGTGTACAGGCTCATCCGCG
ACTCCCGCTACTCAGAGCGGTGCAGCTACTGAGCGCAGAGCTGCAGGGAGTCCCGCAGCCGCGCAGG
GCTGTCGCTGCTGGCCTACTGCTACTACCGCTGCAGGAGTTCGAGCTCGCTGCGGAGTGTATGAGCAA
CTGAGCCAGATGCATCCGGAACCGAGCAGTACCGCTGTATCAGGCCAGGCGCTCTACAAGGCCTGCC
TGTATCCAGAGGCCACCCGGTGCCTTCTCCTGGACAACCCACCTATCAGACTCGAGTCTTCTGCT
GCAGGCTGCTATCAAGTACAGTGAAGGGGATCTCCAGGAGCCAGGAGCTTGGTGGAGCAGCTGCTGAGT
GGGAAGCTGGAGAAGACAGTGGAGGGGAGAATGACCCAGATGGTGGTCAACATGGGTGTCTGCTCT
ACAAGGAGGGACACTACGAAGCTGCCTGTTCAAGTCTTAGCGGCCCTGCAGGCTTCTGGTTACCAACC
TGACCTTCTCAACTGGCTTTAGCCTATTACAGCAGTCGGCAGTATGCGCTGCTCTGAAGCATATC
GCTGACATCATTGAGCGTGGCATCCGTCAACACCCAGAAGTGGTGGGATGACCATGAAGGCATTG
ATATTCGTAGTGTCCGCAACACTGCAGTCTTCATCAGACTGCTCTGATCGAAGCCTTCAACCTCAAGGC
AGCCATAGAGTACCAGCTGAGAACTTTGAGGTAGCCCAAGAAACCTCACTGACATGCCACCCAGAGCA
GAGGAAGAGTTGGACCCCGTGACCTGCACAACAGGCTCTGTGAACATGGATGCCAAGCCACAGAGG
GCTTTGAAAAGCTCCAGTTTCTGCTCCAGCAAAACCCCTCCCCCAGAGACCTTTGGCAACCTGCTGCT
GCTTTACTGTAATATGAGTATTTGACCTGGCAGCTGATGTTCTAGCAGAAAACGCCACTTGACTTAC
AAGTTCCTCACACCCTATCTCTATGACTTCTGGATGCCATGATCACTTGCCAGACAGCTCTGAAGAGG
CTTTCATTAACCTTGATGGGCTAGCAGGTATGCTGACTGACCAGCTCAGGAGACTCACTAAACAAGTTCA
AGAAGCAAGGCATAACAGAGACGATGAAATTATTAAGAAAGCTGTGAATGAATATGACGAAACCTTGAG
AAGTATATCCCTGTACTGATGGCCAGGCAAAATCTACTGGAATCTTGAGAATTATCAAAATGGTGGAGA
AGATCTCCGCAATCTGTGGAATTCTGTAATGACCATGATGTGTGGAAGCTGAATGTGGCCATGTCTCT
CTTCATGCAGGAAAACAAGTACAAGGAGGCCATTGGTTTCTATGAGCCATTGTCAAGAAGAATTATGAT
AACATCTGGGTGTCAGTGCTATCGTGTAGCTAACCTCTGTGTTTCTACATTATGACAAGTCAAAATG
AGGAAGCTGAGGAGCTGATGAGGAAGATTGAGAAAGAGGAGGAACAACCTCTCTATGAAGACCCGGACAA
GAAAATCTACCACCTTTGCATTGTGAATTTGGTATAGGAACCTTTATGTGCCAAAGGAAACTATGAC
TTTGGCATCTCTCGGTTATCAAGAGCCTGGAGCCTTACCATAAGAAGCTGGGAAGTATACGTGGTATT
ATGCCAAAAGATGCTTCTGTCCTTACTGGAGAACATGTCAAAGCACATGATAGTTCTTTGTGATGGTGT
TGTTCAAGAATGTGTCCAGTTTCTAGAGCACTGCGAGATTTTTGGTAGAAACATCCCTGCTGTTCTCGAA
CAACCCTTGGAGGAAGAAAGAAATTCACACTGGAAGAATACAGTCACTTATGAGTCCAGACTGTTGAAAG
CTTTGATTTATGAGATCATAGGGTGAATATG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR217077 representing NM_001271494
Red=Cloning site Green=Tags(s)

```
MAWQSSSKVPDGEFTAVVYRLIRDSRYSEAVQLLSAELQGSRSRAGLSLLAYCYRRLQEFELAAECYEQ
LSQMHPLEQYRLYQAQALYKACLYPEATRVAFLLDNPTYQTRVLRQAQAIKYEGLPGARSLVEQLLS
GEAGEDSGGENDPDGLVNMGCLLYKEGHYEAACSKFLAALQASGYQPDLSYNLALAYSSRQYAPALKHI
ADIIERGIRQHPELGVGMTIEGIDIRSVGNTAVLHQTALIEAFNLKAAIEYQLRNFEVAQETLTDMPPRA
EEELDPVTLHNQALMNMDAKPTEGFQKLLQQLNPFPPETFGNLLLLYCKYEFDLAADVLAENAHLY
KFLTPYLYDFLDAMITCQTAPEEAFIKLDGLAGMLTDQLRRLTKQVQEARHNRDDEI IKKAVNEYDETLE
KYIPVLMQAQAIYWNLENYQMVEKIFRKSVEFCNDHDVWKLNVAVHVMQENKYKEAIGFYEPVVKKNYD
NILGVSAIVLANLVCVSYIMTSQNEEAELMRKIEKEEQLSYEDPKKIYHLCIVNLVIGTLYCAKGNVD
FGISRVIKSLEPYHKKLGTDWTWYAKRCFLSLEENMSKHMIVLDCDGVVQECVQFLEHCEIFGRNIPAVLE
QPLEEERIHTGKNTVTYESRLLKALIYEIIGWNM
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001271494

ORF Size: 1992 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001271494.1](#), [NP_001258423.1](#)

RefSeq Size: 3455 bp


RefSeq ORF: 1995 bp

Locus ID: 100362431

Cytogenetics: 3q23

MW: 76.1 kDa

Product images:



Circular map for RR217077