

## Product datasheet for **RR217083**

### **Pde4c (NM\_001276765) Rat Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Pde4c (NM_001276765) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Pde4c
Synonyms:	Dpde1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide  
Sequence:

>RR217083 representing NM\_001276765  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCCGGATCGCC

ATGAGCGGAGTGGCAGCACCTGTCTTCTGTGACTGAGCGAGTCCGAGAGCCAGTGGATTCTGGAG  
TGGCCCTGTGAGCCACCTAGGTGGGGTGTGATCCTGAGACGGTTCTCAGGGACCCTGCTGCTGCCTCC  
GCTGTCCAGCCGACTCGGGTCTCTGGGAGGCGGAGTCTGCTGCCCCCGTGGTTTTACCATCGGATCC  
CGAGGGACAGAGAGAAGCCATGGCTCAGCTCAGAGCAGCTTTGACCTGAAAAATGGGCTCCCGAGCGCA  
GGGGTCCCCTGGATGCTCAGTCCGGTTCGAGTCTCGGCCGCGCCCTCAAACCCCTGTGCATCACGGCCA  
GCGGCGGAGTCTTCTGTACCGCTCGGACAGTGACCATGAGCCTTACCCAAGGCTGTGTCCCGGACC  
TCCTCGGCGGCCAGCGACCTACATGGAGAAGACATGATTGTCACACCCTTCGCCAGGCTCTGGCCAGTC  
TCAGGACGGTCCGTAGCAACGTGGCTGCCCTTGACATGGCGCGGCAGTGCCACTAGGCAGGCGCTCT  
GGGACCCTCCGAGAGCAGCCAACAGGCCGCCAGCAGAGGAATCTGGGCTGCAGCTGGCGCAGGAG  
ACACTGGAGGAGCTGGACTGGTGCCTGGAGCAGCTGGAGACACTGCAGACTCGTGCCTCGGTGGGCGAGA  
TGGCTTCAACAAGTTCAAGCGCATGCTCAACCGTGAGTTGACCCACCTGTCCGAAAACGACCGCTCAGG  
GAACCAAGTGTGAGATACATCTCGCAGACTTCTCGACAGCAGGCTGAGGTGGAGCTGCCGGCTCCG  
CCCACGGAAGACCACCATGGCCATGGCCAGATCACTGGGCTGCGGAAGTCTTGCCACACCAGCCTTC  
CTACAGCTGCCATCCCGCGCTTTGGGGTCCAGACAGACCAGGAGGAGCAGCTGGCTAAGGAGCTTGAAGA  
CACCAACAAATGGGGATTCGATGTGTTAAAGTGGCTGAGCTGAGTGGGAACCGGCCCTCACTGCTGTT  
ATCTTCAGCGTTCTCCAGAACGGGACCTGCTCAAGACGTTCCAGATCCCGGCGGACACGCTGCTGAGGT  
ACCTGCTGACCTCGAAGGGCACTACCACTCAACGTTGGGTACCACAACAGCATACATGCTGCCGACGT  
GGTGCAATCCGCGCACGTGTTGCTGGGCACCCAGCCCTGGAGGCTGTGTTACGGACCTAGAAGTCTTG  
GCCCCATCTTTCGCTGCGCTATACACGATGTAGACCACCCAGGTGTCTCCAACCAAGTTTCTCATCAACA  
CCAACCTCGGAGCTGGCACTGATGTACAATGACTCATCGGTGCTAGAGAACCATCACCTCGCGGTAGGCTT  
CAAGCTCTGCAGGGGAAAACCTGCGACATCTTCCAGAACCTCAGCACCAAGCAGAAAACCTGAGTCTGCGC  
AGGATGGTCATAGACATGGTACTGGCCACAGACATGTCCAAGCACATGAGCCTCTGGCTGACCTCAAGA  
CCATGGTGGAGACCAAGAAAGTACTAGCCTTGGCGTCTGCTCTTGACAACCTACTCTGACCGCATCCA  
GGTCTCCAGAGCCTGGTGCCTGCGCCGACCTCAGCAACCTGCCAAGCCACTACCCCTCTACCGCCAG  
TGGACGGAGCGCATCATGGCTGAGTTCTTCCAGCAGGGTGACCGGAACGTGAGTGGGCTTGGACATCA  
GCCCCATGTGCGACAAGCACACGCTCGGTGGAGAAATCCAGGTGGGATTCATTGACTACATCGCTCA  
CCCATTGTGGAAACTTGGGCCGACCTGGTGCACCCCGATGCCAGGAGCTGCTGGATACCTTGGAGAC  
AACAGAGAGTGGTATCAGAGTAGGGTCCCTGCAGCCCTCCACACGCCATTGGCCCTGACAGGTTCAAGT  
TTGAGCTGACCCTGGAGGAGACAGAGGAAGAGGAGGAAGAGGATGAAAGGCAC

ACGGTACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RR217083 representing NM\_001276765  
Red=Cloning site Green=Tags(s)


```
MRRSGTTL S FLL TERVREP V DSGVAPVSHLGGGVILRRFSGTLLL P LSSRLGSSGEAESAAPVVF TIGS
RGT ERSHGSAQSSFDLENGLP SGRGPLDAQSGSSLGRALQPPVHHGQRRESFLYRSDSDHEPSPKAVSRT
SSAASDLHGEDI V T PFAQVLA SLRTVRSNVAALAHGAGS ATRQALLGTPPQSSQAAPAEESGLQLAQE
TLEELDWCLEQL E T L Q T R R S V G E M A S N K F K R M L N R E L T H L S E T S R S G N Q V S E Y I S Q T F L D Q Q A E V E L P A P
P T E D H P W P M A Q I T G L R K S C H T S L P T A A I P R F G V Q T D Q E E Q L A K E L E D T N K W G F D V F K V A E L S G N R P L T A V
I F S V L Q E R D L L K T F Q I P A D T L L R Y L L T L E G H Y H S N V A Y H N S I H A A D V V Q S A H V L L G T P A L E A V F T D L E V L
A A I F A C A I H D V D H P G V S N Q F L I N T N S E L A L M Y N D S S V L E N H H L A V G F K L L Q G E N C D I F Q N L S T K Q K L S L R
R M V I D M V L A T D M S K H M S L L A D L K T M V E T K K V T S L G V L L L D N Y S D R I Q V L Q S L V H C A D L S N P A K P L P L Y R Q
W T E R I M A E F F Q Q G D R E R E S G L D I S P M C D K H T A S V E K S Q V G F I D Y I A H P L W E T W A D L V H P D A Q E L L D T L E D
N R E W Y Q S R V P C S P P H A I G P D R F K F E L T L E E T E E E E E E D E R H
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_001276765

**ORF Size:** 2013 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001276765.1](#), [NP\\_001263694.1](#)

**RefSeq Size:** 2161 bp

**RefSeq ORF:** 2016 bp


**Locus ID:** 290646

**Cytogenetics:** 16p14

**MW:** 74.2 kDa

**Gene Summary:** putative nucleotide phosphodiesterase [RGD, Feb 2006]

### Product images:



Circular map for RR217083