

## Product datasheet for **RR217085**

### Srp72 (NM\_001170601) Rat Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Srp72 (NM_001170601) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Srp72
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide  
Sequence:

>RR217085 representing NM\_001170601  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCGAGCGGCGCAGCGGTGGGTGTCACTCCCTGCGCTGTGGAGTGAGGTGAACCGTTATGGGCAGA  
ACAGCGACTTCACGCGTGCCCTCAAGACCGTCAACAAAATATTGCAGATCAACAAGGATGATGTCACGGC  
CTTGCACTGTAAAGTTGTTTGCCTTATCCAGAATGGAAGTTTCAAAGAAGCCTTGAATGTCATCAATACT  
CACACCAAAGTGTAGCCAATAACTCTCTGTCTTTGAGAAGGCATATTGCGAGTACAGGTTGAACAGAA  
TCGAGAATGCCTTGAAGACAATAGAAAGTGCCAACCAGCAGACAGACAAAAGTGAAGGAGCTTTATGGACA  
AGTGCTGTACCGCTGGAACGCTACGATGAGTGCCTGGCCGTGTACCGAGATCTTGTCCGAAACTCCCAG  
GATGACTACGATGAGGAGAGGAAGACAAACCTGTCAGCAGTCGTTGCCGCTCAAAGCAACTGGGAAAAAG  
TGTTCCCGAGAACTGGGTCTCCAAGAAGGCACACGAACTCTGTTACAACGCTGCGTGTGCACTGAT  
AGGGCAAGGCCAGCTGACCCAGGCAATGAAAATCTACAAAAGGCTGAAGATCTTGTGCGCGTTTCAATT  
TCAGAAGATTCAGACGGGACTGAGGAAGACCCTCAGGCCGAGCTGGCCATCATTACCGTCAAATGGCTT  
ACATCATGCACTTCAGGGGCGCATCGAGGAGGCTCTGCACTTTACAATCAGATAATAAAAAGTGAAGCC  
AACAGATGTGGCACTTCTTGCTGTAATCGCAAATAACATCATCACCATTAATAAGGACCAAAAATGTCTTT  
GACTCCAAGAAGAAAGTGAATTAACCAATGCAGAAGGAGTAGAGTTAAGCTCTCCAAGAGACAGCTTC  
AAGCCATAGAGTTCAACAAGGCTTTGCTTGCCATGTATACAAACCAGGCAGAACAGTGCCGAAAAATAGC  
TTCTAGCTTGCACTCCAGAGTCTGAGCATCTCTGCCTGTGCTGATCCAAGCTGCCAGCTCTGTGCTG  
GAAAAGCAGCATACAAAGGCAATAGAGCTGCTCAGGAATTTCTGATCAACATCCAGAAAATGCAGCTG  
AAATCAAGCTTACCATGGCTCAATTGAAAATTTCCCAAGTAATATTTCCAAAGCGTGCCTAATACTGAG  
AAGCATAGAGGAGTTAAAGCATAAGCCAGGCATGGTATCCGCGTTGGTGACCATGTACAGCCATGAGGAA  
GACATTGATAGTGCCATTGAAGTCTTACACAGGCTATCCAGTGGTATCAAAGCCATCAGCCCAAGTCCC  
CTGCTCATTTGCTTGTATAAGAGAGGCCGCAAACCTTCAAGCTCAAATATGGACGGAAGAAGGAGGCAGT  
TAGTGACCTAGAGCAGCTGTGGAACAAAATTCAAAAGATATTCACACGCTAGCCAGCTTATCTCTGCT  
TACTCACTTGTGGATCCAGAGAAAGCAAAGGCTCTTAGTAAACACCTGCCCTCATCGGATAGTATGTCTC  
TAAAGGTAGATGTTGAGGCTCTGAGAATTCCTCTGGTGTACGTACATTCGGAAGAAGGGTGGAAAAGT  
TACTGGAGATAATCAACCAAAGGAACAAGGACAGGAGATTTGAAGAAGAAGAAGAAAAAGAAAGGA  
AAACTGCCTAAGAATTATGACCCAAAAGTTACCCAGATCCAGAAAGATGGCTACCAATGCGAGAGCGCT  
CTTACTACCGAGGAAGAAAGAAGGGCAAAAAGAAGGACCAGATCGGAAAAGGAACACAGGAGCGACCCG  
AGGAGCCTCCTCTGAGCTGGACGCCAGTAAAAGTGTGAGCAGCCACCCACCTCCCAAGACCTGGCAGT  
GCGGCAACGATGGCTCCTCTACAAGTAACATTGTACCCCAAGACACCAGAAACCCGACGGGGCCCGAG  
CAACGAAAAAGAAGCAGCAACAGAAGAAGAAGAAGGAGGGAAAAAGTGGTTGG

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC  
TGGATTACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RR217085 representing NM\_001170601  
 Red=Cloning site Green=Tags(s)


MASGGSGGVSPALWSEVNRYGQNSDFTRALKTVNKILQINKDDVTALHCKVVCLIQNGSFKEALNVINT  
 HTKVLANNLSFEKAYCEYRLNRIENALKTIESANQOTDKLKYQVLYRLERYDECLAVYRDLVRNSQ  
 DDYDEERKTNL SAVVAAQSNWEKVVPENLGLQEGTHELCYNAACALIGQGQLTQAMKILQKAEDLCRRSF  
 SEDSDGTEEDPQAE LAI IHQMAYIMQLQGRIIEALQLYNQI IKLKPTDVAL LAVIANNIITINKDQNVF  
 DSKKKVKLTNAEGVEFKLSKRQLQAIIEFNKALLAMYTNAEQCRKIASLSQSQSPEHLLPVL IQAAQLCR  
 EKQHTKAIELLQEFSDQHPENAAEIKLTMAQLKISQGNISKA CLILRSIEELKHKPGMVSALVTMY SHEE  
 DIDS AIEVFTQAIQWYQSHQPKSPAHL SLIREAANFKLYGRKKEAVSDLEQLWKQNSKDIHTLAQLISA  
 YSLVDPEKAKALSKHLPSSDSMSLKVDVEALENSPGATYIRKKGKVTGDNQPK EQGQDLKKKKKKKKG  
 KLPKNYDPKVTDPERWLP MRERSYYRGRKKGKKKDKQIGKGTQGATAGASSELDASKTVSSPPTS PRPGS  
 AATMASSTSNIVPPRHQKPAGAPATKKKQKQKQKGGKSGW

SGP TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

**Restriction Sites:**

SgfI-RsrII

**Cloning Scheme:**



**ACCN:** NM\_001170601

**ORF Size:** 2013 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** NM\_001170601.1, NP\_001164072.1

**RefSeq Size:** 3597 bp


**RefSeq ORF:** 2016 bp

**Locus ID:** 498351

**Cytogenetics:** 14p11

**MW:** 74.7 kDa

**Product images:**



Circular map for RR217085