

Product datasheet for **RR217097**

Lrrc43 (NM_001170396) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Lrrc43 (NM_001170396) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Lrrc43
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RR217097 representing NM_001170396
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGGCGACCAGAGGGTCCAGCGCTAGCGACGAAAGGCATCTCGGGCCGAGGACGTACCAGGGACCTTGA
 GCACGGCTGTGCGTGAACACCTGCGGAAGCTGTGTCTGCGCGAATTCCTGTTGGCGTCGGAAGCTGGAA
 TAAATCACGCTTTCTTCCTCGACAATATCGAACCTGGAGGGAATTGATCCCCAAGGAGGAGAGACAGTG
 GGTCCCAGGAGGAGACCGTGGAGGCCCTGCTGGGCTGGTCCGCAGCAGCCACTCTCCATGGGCTATGC
 TGAAGGATTCCAGCGCAGAGGACCGTTTCTGAGAGAAGTGGCCATCCAGAATCCACTGACGATCAAAGA
 CACTTTCTTACTCTACTTCCGGTGCCTGCGGGTGGTGGACAAGGGGGTGAAGTCTGGTGGATAAAGAA
 ATCTTGAAATTTCTAAAGCTCGAGGAGTTGGTCTGAGTGCCAATAAAATCAAGGACATCGATGCCAACA
 ATCTGCCCCCAACTGAAGGTGCTGGAGCTGTACGGCAATTTGATCGCCAGCATGGAGTGTCTGTGCTC
 ACCCCACCCCAAACTGCAGCACCTGGGGCTAGGCCACAACAAGTTGCTGGGCCCTTTGGAAAGTCTG
 TATGTCACCTCCAACAACCTGGCCTCTTCTCGTCTCCCTGGACCTGGGCTTCAATGACTTGACAGACCTGC
 AGAGCATGATACTGGGTCTCAGTACCTTAAAGTGCCTGAGACTCCTTGTGCTGCAGGGGAACCCACTGGC
 TCTGGTGCCTTACTACCGTGGCTTTACCGTGGACAGCCTGGCGCGCCTCTGTGTGCTGGACGACATCACT
 GTGTCTCCAGCGAGAAGCACCAGTTCGGGGACTCAACAGTCATGGCGACCTGTTGGCAGTGGAGGCC
 AGTTTGTGGTGACCATTGGAATGTGAGGGGGTGGGACAGCTCTATCTTGGACCCAGAGCCAGGACC
 CGATGGCCCTTTCATCAGCTACAGCTATTACGTGACCTATGACTTTGTGGAAGATGAGGAGAAGGAAGGG
 ATTACAAGTGGGATGGAGGAGTTACACCCTGACAGCGTTCGGATGAGATCGACAAGCCCTCCTCTGGCA
 TGGACGAGGAGGAGGAGGAGGAGGAAGAGGAGGAAGACCCAGAGAGGTTCTCAGTGTTCAGAG
 CCAAGCCACCACCGCCCGTCCGCCAGGGTGGCCGGCACAGGTTTCGTCAGAGTCTCCTGAGGTGATG
 TCCAGGGAGCTGTCTGAGTATATAGTTAAGGAGATGAATCAGATGGCCGAGGACTCCGGAGAATCTGGGA
 TCACAGAAGTAGATGAGTCAGAGACAAGCATTTCCATCCACTCGGCCCCCTGCCACAGTCCATTGACTC
 TTCCGAGGAGCTGGCGAAGCTGCGGCCACGGCTAGACGTTGCGCTCTGTCCATCCCAGGGACAGTGCTC
 TTCAGCACCATGCACAAGCCCTGGGCAGATGTCATCCCTTGACCTATGAGATGAAGCACACCCTAAAGG
 AGCTGATAAGGGTCAAGGCCTTCTGTGTCTGGGACCACTGTGAGCATCGTAGAGGAGAAAATCATCTC
 CTGGCCTGTGGTACCAACTCCTGTTGAGAGTCCCTTGCCTGCCAAGAAGGGGAAAGGCGACAAGAAAAAG
 AAAGAGGAGCCTGCTAAGGACAAAGGCCGAGGAAGAAAAGGAGCCACCAGAGAGCTGCGCCAGGATC
 CACCTGTGCTGAGGGTGTGGGCAGCGGACAGGTGTACCTGGAGCCCTGCTGGCAGGGGAGCCAGTGGT
 GTCCACTGTGTGAATTTGGGGTGGCCCGCACCTGGAGACTGACAGGCTCACATATGCCAGGGATTCA
 AAGAAGGTAAGAAAGGCTCAAAAAAGATAAGTCAAAAGCGCGGCTCCAACCGTAGAAAGCGGGTACC
 AGCCAGAGCCCTGTGGTGGAGGTCCAGATCCAGTTGCACCAGTACCGCTCGGTGGAGGAAGCGTTCTCT
 TAGCTTGATTGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR217097 representing NM_001170396
Red=Cloning site Green=Tags(s)


MATRGSSASDERHLGPEDVPGTLSTAVREHLRKLCLREFPCGVGSWNKSRFLPRQYRTWRELIPKEEETV
 GPEEETVEALLGLVRSSHSPWAMLDSSAEDRFLRELAIQNPLTIKDTFFYSYFRCLRVVDKGVSLVDKE
 ILKFLKLEELVLSANKIKDIDANNLPPTLKVLELYGNLIASMECLCSPPPPNLQHLGLGHNKLLGPLESL
 YVTSNNWPLLVSLLDLGFNDLTDLQSMILGLSTLKCLRLVLQGNPLALVPYYRGFTVDSLARLCLVDDIT
 VSPSEKHQFRGLNSHGDLLAREAQFVVTIGNVRGVLDSIILDPEPGDPGPFISYSYVYTYDFVEDEEKEG
 ITSGMEELHPDSVLEIDKPSSGMDEEEEEEEEEEDRQEVLSVFQSQAHHRPVRQGRHRFRPESPEVM
 SRELSEYIVKEMNQMAEDSGESGITEVDESETSISIHSAPLPQSIDSSEELAKLRPRLDVRLCPSPGTVL
 FSTMHKPWADVIPCTYEMKHTLKEKIRVKAFLLSGTTVSIVEEKIISWPVVPTPVESPLAKKGGKDKKK
 KEEPAKDKGPRKKKEPPRELQDPPVLRVLGSGQVYLEPLLAGEPVVSTVCNFGVARTLETDRLYARDS
 KKVKKGSKDKSKAAAPTVESGYQPEPLSVEVQIQLHQYRSVEEAFSLSID

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001170396

ORF Size: 2043 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001170396.2](#), [NP_001163867.1](#)

RefSeq Size: 2097 bp


RefSeq ORF: 2046 bp

Locus ID: 288751

Cytogenetics: 12q16

MW: 76.3 kDa

Product images:



Circular map for RR217097