

Product datasheet for **RR217100**

Ablim3 (NM_001191698) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ablim3 (NM_001191698) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ablim3
Synonyms:	RGD1565118
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RR217100 representing NM_001191698
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAACACGAGCATTCCCTTACCAGCAGAGTCCGTACAGTCCCCGGGGCAGCTCCAACGTCATCCAGTGCT
 ACCGCTGCGGAGACACCTGCAAAGGAGAGTGGTCCGTGTCCACAACAACCACTTCCACATCAGATGCTT
 CACTTGTCAAGTATGTGGATGTGGCCTGGCCAGTCGGGCTTCTTCTCAAGAACCAGGAGTACATCTGC
 ACCCAGGATTACCAGCAGCTCTATGGCACCCGCTGTGACAGCTGCCGGGACTTCATCACGGGGGAGGTCA
 TCTCCGCCCTGGGCCGCACCTACCACCCTAAGTGTCTCGTGTGCAGCCTGTGCAGGAAGCCTTCCCTAT
 TGGAGATAAGGTGACCTTCAGTGGTAAAGAATGTGTCTGTGACAGCTGCTCCAGTCGATGACCAGCAGC
 AAGCCAATCAAGATCCGTGGACCAAGCCACTGTGCTGGTGCAAAGAGGAGATTAACACGGCCAGTCAC
 TCCTGGCACTGGACAAGCAATGGCACGTCAGTGTTCAAATGCCAGACCTGTAGCGTCATCCTCACTGG
 GGAATACATTAGCAAGGATGGTGTCCGTACTGCGAGTCTGATTACCACGCCAGTTTGGCATCAAATGT
 GAGACTTGTGACCGGTACATCAGTGGCAGGGTCTTGGAGGCAGGAGGAAAACACTACCACCTACCTGTG
 CCAGATGTGTGGCCTGCCACCAGATGTTCACTGAAGGGGAGGAGATGTACCTCACAGTTCTGAGTTTTG
 GCACCCGATCTGCAAGCAGGCAGCCAGGGCAGAGAAGAAGCTCAAGCATAGACGGACATCTGAAACTTCC
 ATCTCACCCCCGGATCCAGCATCGGGTGCCTAACCGAGTCACTGTGCTAAAGTGATAATGAGATCC
 TTAATTACAAAGACCTGGCGGCTCTCCCAAGGTTAAGTCCATCTACGAGGTGCAACGTCCCAACCTCAT
 TTCTACGAGCCTCATTCCAGATACACCTCCGACGAGATGCTGGAGAGATGCGGTTACGGAGAGTCACTG
 GGAACATTATCTCTTACTCCCAGGACATCTATGAGAACTGGACCTCCGGCAGAGACGGCCCTCCAGCC
 CGGGATACATCGATTCCCCCCTACAGTCAGCAAGGCATGTCTCCACCTTCTCCGCTCACCTCATCA
 TTAATACCGCTCAGGGCCTGAGAGTGGCCGGAGCTCTCCGTACCATAGCCAGTTAGATGTGAGGTCCTCC
 ACTCCGACCTCTTACCAGGCTCCTAAGCACTTTCATATCCCAGCTGGAGAAAGTAAACATCTACCGAAAAC
 CCCCAGTCTACAAACGGCATGGTATTTGTCCACAGCAACCAAGAGCAAGACGAGTGAAGACATCAGCCA
 GGCCCTCAAGTACAGCCAGCCTACTCACCAGACCCCTACTATGCTTCAGAGTCTGACTATTGGACCTAC
 CATGGGTCCCCAAAGTGGCCAGAGCCAGAAGTTCTCATCTGGAGGAGAGGAAGACGATTTTGACCGCA
 GTATGCACAAGCTCCAGAGTGAATTGGCCGGCTGATTCTGAAGGAAGAAATGAAGGCACGATCAAGCTC
 CTATGCGGATCCTTGGACCCCTCCGCGGAGCTCCACCAGCAGTCGGGAGGCCCTGCACACCACTGGCTAT
 GAGATGTCCTTCAATGGCTCCCCCGTTACACTACCTGGCTGACAGTGACCCTCTCATCTCCAAGTCTG
 CCTCACTGCCTGCTACAGGAGAAATGGGCTGCACAGGACACCCAGTGGCGACCTTCCACTACGACAG
 CATGAACGCAGTCAACTGGGCATGCGAGAGTACAAGATCTACCCATATGAGCTGCTGCTGGTGACCACA
 AGAGGAAGGAACCGGCTGCCAAAGGATGTGGACCGAACTCGTTAGAGCGCCACCTGTCCCAGGAGGAGT
 TCTACCAAGTCTTCGGCATGACCATCTCAGAGTTCGAGAGGCTAGCCCTCTGGAAGAGGAATGAAGTAA
 GAAGCAAGCCCGGCTGTTC

ACGCGTACGCGGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR217100 representing NM_001191698
Red=Cloning site Green=Tags(s)


MNTSIPYQQSPYSPRGSSNVIQC YRCGDTCKGEVVRVHNNHFHRCFTCQVCGCLAQSGFFFKNQEYIC
 TDQYQQLYGRCDSCRDFITGEVISALGRTYHPKCFVCSLCKRPFPIGDKVTFSGKECVCQTCSQSMSS
 KPIKIRGPHSHCAGCKEEIKHGQSLLDKQWHVSCFKCQTCVILTGEYISKDGVYPYCESDYHAQFGIKC
 ETCDRYISGRVLEAGGKHYHPTCARCVRCHQMFTEGEEMYL TGSEVWHPICKQAARAEEKLKHRRTSETS
 ISPPGSSIGSPNRVICA KVDNEILNYKDLAALPKVKSIYEVQRPNLISYEPHSRYTSDEMLERCGYGESL
 GTLSPYSQDIYENLDLRQRASSPGYIDSPTYSQQGMSPTFSRSPHHYRSGPESGRSSPYHSQLDVRSS
 TPTSYQAPKHFHIPAGESNIYRKPPYIKRHGDLSTATKSKTSEDISQASKYSPAYSPDPYASESDYWTY
 HGSPKVPARRFSSGGEEDDFDRSMHLQSGIGRLILKEEMKARSSYADPWTPPRSSTSSREALHTTGY
 EMSFNGPSRSHYLADSDPLISKSASLPAYRRNLHRTPSADLFHYDSMNAVNWGMREYKIYPYELLLVTT
 RGRNRLPKDVRTRLERHLSQEEFYQVFGMTISEFERLALWKRNELKKQARLF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001191698

ORF Size: 2049 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001191698.1](#), [NP_001178627.1](#)

RefSeq Size: 2833 bp


RefSeq ORF: 2052 bp

Locus ID: 307395

Cytogenetics: 18q12.1

MW: 77.7 kDa

Product images:



Circular map for RR217100